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HSW HENKE-DENT 300 A Trio HSW HENKE-DENT 300 E Trio

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Stand: 09.09.2019


HSW HENKE-DENT 300 A Trio HSW HENKE-DENT 300 E Trio

Intended use: Dental cartridge syringe for infiltration- and block anaesthesia.

This syringe with aspiration function is used together with a sterile dental needle and a cartridge, prefilled with an anaesthetic, primarily for local anaesthesia especially in dentistry (both are not included in scope of delivery).

Assembly of: HSW HENKE-DENT 300 A Trio (chromed brass)
HSW HENKE-DENT 300 E Trio (stainless steel)

The HSW HENKE-DENT 300 A / E Trio is a non-detachable injection syringe.

Cleaning & Sterilisation

Prior to the initial use of the product and for each further use the product must be cleaned and sterilized (Syringes will be delivered clean and non-sterilized):

Preparations at the place of installation:	Remove obvious dirt particles from medical devices immediately after utilization. It is not allowed to use fused agents or hot water (> 40°C) for cleaning as these agents lock residues and impede an effective purification process.
Transport:	Secure storage as well as secure transport in closed boxes is recommended in order to avoid any damage to the medical device as well as environmental contamination.
Manual pre-cleaning:	The medical device must be cleaned by means of a soft multi-purpose brush and cold tap water until all obvious residues and defilements are removed. Lumen, drill holes and thread turns must be flushed for at least 10 seconds at a pressure of 3.8 bars by means of a water pistol.
Cleaning	Medical devices must be put in a screen basin on the plug-in mode. Then start the cleaning process. <ol style="list-style-type: none"> 1. 4 minutes pre-washing by means of cold water 2. Draining 3. 10 minutes pre-washing with 0,5% Neodisher Mediclean, Dr. Weigert (Hamburg) 4. Draining 5. 6 minutes neutralization with 0,1 % Neodisher Z, Dr. Weigert (Hamburg) (> 40°C) 6. Draining 7. 3 minutes flushing with VE-water (< 40°C) 8. Draining
Desinfection:	The automatic thermal disinfection must be carried out according to the national requirements regarding the A0-coefficient (see ISO 15883).

Drying:	Drying of the exterior parts of medical devices by means of the drying cycle of the cleaning / disinfection equipment. Additional manual drying is possible by means of a lint-free cloth. Hollow space of medical devices must be dried by means of sterile compressed air.
Functional test, maintenance	Afterwards, a visual appraisal regarding purity is necessary. The maintenance and functional tests must be carried out according to the instruction manual. If necessary, the reprocessing procedure must be re-run until the medical device is obviously clean.
Packaging	Packaging of medical devices for sterilization must comply with standards according to ISO 11607 and EN 868.
Sterilization (autoclaving):	Sterilization of products by means of the fractionated pre - vacuum method (according to ISO 13060 / ISO 17665) in compliance with the individual national requirements. 5. Fractionated pre - vacuum method (3-fold) 6. Sterilization temperature of 134°C 7. Minimum exposure time: 3 minutes (full cycle) 8. Drying time: at least 10 minutes
Storage:	Storage of sterilized medical devices must be carried out in a dry, clean and dust-free environment at a moderate temperature of 5°C - 40°C.
Information on validation of processing:	The following materials and machines have been used for validation: Cleaning supplies: Neodisher Mediclean (alkaline); Dr. Weigert; Hamburg Neutralization supplies: Neodisher Z; Dr. Weigert; Hamburg Cleaning / disinfection equipment: Miele G 7836 CD
Additional instructions:	According to the Ordinance on Installation, Operation and Use of Medical Devices (MPBetreibV in German), the operator is responsible for validating any reprocessing procedures, even if the above described chemicals and machines are not available. The operator must ensure that any reprocessing procedures – including resources, material and staff - are capable of achieving the required results. National legislation as well as the state of technology requires compliance with validated procedures.

Application

Infiltration anaesthesia:

In infiltration anaesthesia, the local anaesthetic is infiltrated into the tissue of the target area to be anaesthetised.

Block anaesthesia:

In the lower jaw, mainly block anaesthesia is used. The local anaesthetic is inserted as near as possible to the nerve trunk and the entire supply area of these nerves is anaesthetised.

Contra indication:

The contra indications of the anaesthetic used in each case apply.

Details

The HSE-DENT 300 A/E Trio can be used for all standardized glass-cylinder cartridges 1,7 ml or 1,8 ml with hole plugs.

When inserting or removing the cylinder cartridges, the aspiration hooks must be completely retracted (left-hand rotation of the knurled screw) - important to prevent damages to the aspiration hooks.

When angulating of the rear part of the syringe, as well as during the return, the piston rod must be completely retracted during the entire angulating procedure. Otherwise, there is a risk of damaging the piston rod or the aspirating hooks.

Do not unfix or remove the existing screws (angulating mechanism 2 x and fingerrest 1 x) before reprocessing.

In order to ensure a long lifetime of the product, we recommend the strictly observe of our instructions for use, as well as careful handling and preparation.

Preparation of device

8. Pull the plunger rod back to the stop by overcoming the spring resistance and angulate the rear part of the syringe.
9. Insert a suitable cartridge (1,7 ml or 1,8 ml) with the appropriate anaesthetic from behind into the cartridge holder. Make sure that the head of the cartridge is facing forward during the procedure of inserting the cartridge.
10. To close the syringe, fold the cartridge holder back to the stop by overcoming the spring resistance and angulate the rear part back again.
11. Now press the closed aspiration hooks up to the stop into the inserted cartridge and solidly connect the piston rod with the hole plug. The connection is made by right hand rotation of the knurled screw (=extending of the hooks). As soon as the hooks are swung out in the end position, the locking arrest secures the connection and aspiration is possible.
12. Insert a sterile dental injection needle, according to the requirements, into the head of the cartridge holder. Fix the injection needle by screwing it onto the thread of the cartridge holder. Please note that only injection needles with a metric thread will fit on these syringes.
13. The possibility of injection including the possibility of aspiration is constructed now.
If a correct aspiration possibility exists, the function should be tested by slightly pulling back the piston rod.

To check whether the injection needle is correctly mounted, the flow of the anaesthetic should be checked by a slight forward movement of the plunger.

If there is a smooth flow of the anesthetic, the dental cartridge syringe is now ready for application.
14. The removal of the cartridge occurs in reverse order.

Warning notices:

To avoid needle stick, please pay attention to the national and internal specifications and safety instructions.

HENKE-SASS, WOLF GmbH
Keltenstraße 1
78532 Tuttlingen, Germany

www.henkesasswolf.de
info@henkesasswolf.de

Stand: 09.09.2019