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



Instructions for use Dental injection syringes for intraligamental anesthetics (indication)

HSW LIGMAJECT®

HSW HENKE-JECT®

HSW HENKE-JECT® with pressure limit

IMPORTANT: Our syringes are supplied clean but non-sterile and must be cleaned and sterilized prior to each use.

Description:


Dental syringes are used with a sterile dental needle and with a pre-filled anaesthetic cartridge barrel (both not included in delivery).

Warnings:

To avoid needle stick please pay attention to the national and internal specifications and safety instructions. With the HSW LIGMAJECT® respectively the HSW HENKE-JECT® it is not possible to aspirate, they must not be used for infiltration anaesthesia or conduction anaesthesia.

The syringe is delivered with a protective tube. This protective tube is an additional protection if the ampul breaks when there is a higher pressure. Glass splinter in the oral cavity will be avoided.

Function and handling:

- Actuate the brake disc (1.1: HSW LIGMAJECT®) respectively (1.2: HSW HENKE-JECT®) the bolt for releasing the arrested piston rod (2). While you press back the brake disc (1.1) of the HSW LIGMAJECT®, respectively hold the bolt (1.2) of the HSW HENKE-JECT®, you can extract the piston rod (2) to the back. Please note that during the actuation of the brake disc (1.1) respectively of the bolt (1.2) the release lever (4) must be at a passive state. Otherwise a retraction of the piston rod (2) is not possible. After having pressed back the piston rod (2) to the back, release the pressure on the brake disc (1.1) respectively on the bolt (1.2), to prevent a sliding back. Now the piston rod (2) is arrested in the back.
 
- Screw the cartridge holder (3) out of the handle part.
- Then introduce the cylindrical protective tube into the cartridge holder (3).
 
- Now you can introduce the appropriate standard cylindrical cartridge (1.7ml, 1.8ml or 2.2ml) with the corresponding anaesthetic in the unscrewed cartridge holder (3), in which already the protective tube is located. Make sure that during the introduction of the cylindrical cartridge the cartridge head is directed forward.
- Screw the cartridge holder which is now equipped with the protective tube and the cylindrical cartridge, carefully back into the handle part. Please pay attention that the protective tube and the cylindrical cartridge do not slide out during the screwing-in process!
- After that, introduce a sterile dental injection needle corresponding to the requirements into the head of the cartridge holder. Fix the injection needle by screwing it on the thread of the cartridge holder. Note that only needles with a metric thread fit on injection syringes with a metric front part. Injection needles with an inch thread can only be used with injection syringes with an inch.

7. For checking if the injection needle is situated correctly on the cylindrical cartridge, the flow of the anaesthetic should be tested by a slight motion of the release lever (4). If there is a smooth flow of the anaesthetic, the syringe is ready for application. Every motion of the release lever corresponds to 0.2 ml.
8. Please note that the injection must be done slowly.

Instruction for the HSW HENKE-JECT® with pressure limit

Each pressure, which exceeds the data ($90 \text{ N} \pm \frac{7}{9} \text{ N}$ and $120 \text{ N} \pm 12 \text{ N}$) during the application, immediately activates the limiting mechanism whereupon you can continue the injection with a maximum pressure of 90 N respectively 120 N. As soon as the pressure, which is controlled by the release lever, is below the named data, the dentist can control the injection himself.

Sterilization and Cleaning Instructions

Cleaning:

The medical device must be disassembled and cleaned under cold tap water by means of a soft multi-purpose brush until all obvious residues and defilements are removed like blood, saliva and anesthesia solutions.

Disinfection:

The automatic thermal disinfection must be carried out according to national requirements regarding the A_0 -coefficient (see ISO 15883).

Drying:

Drying of the exterior parts of medical devices by means of the drying cycle of the cleaning / disinfection equipment. Additional manual drying is possible by means of a lint-free cloth. Hollow spaces of medical devices must be dried by means of sterile compressed air.

Sterilization (Autoclaving):

Sterilization of products by means of the fractionated pre-vacuum method (according to ISO 13060 / ISO 17665) in compliance with the individual national requirements.

1. Fractionated pre-vacuum method (3-fold)
2. Sterilization temperature of 134°C
3. Minimum exposure time: 3 minutes (full cycle)
4. Drying time: at least 10 minutes

Please note that the protective sleeve (protective tube) lies always offside from the metal components. Reassemble the components only in completely dry condition. Oiling is not necessary. Detailed instructions for automatic preparation see also: www.henkesasswolf.de or can be requested directly from Henke-Sass, Wolf GmbH.

Indication of the method intraligamental anaesthesia (ILA)

Preparations of cavities and crowns.

Endoscopic treatment with the possibility of re-injection in the pulp space.

Single dental retraction in permanent dentition.

Differential diagnosis of unknown pulp affections.

Treatment of patients with instable circulation, because of the small amount of the necessary anaesthetic.

Patients with haemorrhagic diathesis and medicated with Marcumar.

Contraindicated indication of the method intraligamental anaesthesia (ILA)

Anaesthetization of patients with risk after heart attack, cardiac bypass and other ailment.

Treatment of children and handicapped to reduce the risk of postoperative bite injuries.

In the primary dentition possibility of a hypoplasia.

Completion of partial anaesthesia failure concerning conduction anaesthesia.

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