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### HSW Soft-Ject® Dosierradspitze

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## Instruction for use

### HSW Soft-Ject® Micro-dosage wheel syringe


**Intended use:** Micro-dosage wheel syringe for intraligamental injections

The pressure required for the intraligamental injection is built up by moving the dosing wheel with the thumb or fore finger. The dosing wheel transmits the movement to the piston rod and consequently to the gasket of the cartridge. As a result, the anaesthetic is injected.

#### Assembly of the HSW Soft-Ject® dosing wheel syringe

The HSW Soft-Ject® dosing wheel syringe consists of:

- Cartridge holder consisting of a fenestrated sleeve, burst protector and needle connector
- Instrument holder consisting of dosing wheel and toothed piston rod for continuous pressure transmission
- Both parts are connected with a bayonet lock

#### Cleaning & Sterilisation

Prior to the initial use of the product and for each further use the product must be cleaned and sterilised:

<b>Preparations at the place of installation:</b>	Remove obvious dirt particles from medical devices immediately after utilization. It is not allowed to use fused agents or hot water (> 40°C) for cleaning as these agents lock residues and impede an effective purification process.
<b>Transport:</b>	Secure storage as well as secure transport in closed boxes is recommended in order to avoid any damage to the medical device as well as environmental contamination.
<b>Preparations for decontamination:</b>	For further reprocessing, medical devices must be disassembled or opened up as far as possible. For details see separate instructions.
<b>Manual pre-cleaning:</b>	The medical device must be disassembled (see above) and cleaned by means of a soft non-abrasive brush and cold tap water until all obvious residues and defilements are removed. Lumen, drill holes and thread turns must be flushed for at least 10 seconds at a pressure of 3.8 bars by means of a water pistol.

<b>Cleaning:</b>	Disassembled or opened-up medical devices must be put in a screen basin on the plug-in mode. Then start the cleaning process. 1. 4 minutes pre-washing by means of cold water 2. Draining 3. 5 minutes pre-washing at 55° C with 0,5% Neodisher Mediclean, Dr. Weigert (Hamburg) 4. Draining 5. 3 minutes neutralization with warm tap water (> 40°C) 6. Draining 7. 2 minutes intermediate flushing by means of warm tap water (> 40°C) 8. Draining
<b>Disinfection:</b>	The automatic thermal disinfection must be carried out according to the national requirements regarding the A0-coefficient (see ISO 15883).
<b>Drying:</b>	Drying of the exterior parts of medical devices by means of the drying cycle of the cleaning / disinfection equipment. Additional manual drying is possible by means of a lint-free cloth. Hollow space of medical devices must be dried by means of sterile compressed air.
<b>Functional test, maintenance:</b>	Afterwards, a visual appraisal regarding purity is necessary. The assembling of medical devices, maintenance and functional tests must be carried out according to the instruction manual. If necessary, the reprocessing procedure must be re-run until the medical device is obviously clean.
<b>Packaging:</b>	Packaging of medical devices for sterilization must comply with standards according to ISO 11607 and EN 868.
<b>Sterilization (autoclaving):</b>	Sterilization of products by means of the fractionated pre - vacuum method (according to ISO 13060 / ISO 17665) in compliance with the individual national requirements. 1. Fractionated pre - vacuum method (3-fold) 2. Sterilization temperature of 134°C 3. Minimum exposure time: 3 minutes (full cycle) 4. Drying time: at least 10 minutes
<b>Storage:</b>	Storage of sterilized medical devices must be carried out in a dry, clean and dust-free environment at a moderate temperature of 5°C - 40°C.
<b>Information on validation of processing:</b>	The following test instructions, materials and machines have been used for validation:  Cleaning supplies: Neodisher Mediclean (alkaline), Dr. Weigert; Hamburg  Cleaning / disinfection equipment: Miele G 7735 CD with plug-in module, vario-TD-program (without disinfection)  For details see report Cleaning: 17607011411 - 1 Sterilization: 17607010811 - 1
<b>Additional instructions:</b>	According to the Ordinance on Installation, Operation and Use of Medical Devices (MPBetreibV in German), the operator is responsible for validating any reprocessing procedures, even if the above described chemicals and machines are not available. The operator must ensure that any reprocessing procedures – including resources, material and staff - are capable of achieving the required results. National legislation as well as the state of technology requires compliance with validated procedures.

#### Application

Intraligamental anaesthesia allows for a safe and deep short-term analgesia of a single tooth. The nerve ends surrounding the single tooth and pulp nerves are desensitized by injecting the anaesthetics into the periodontal gap.  
The required amount of anaesthetics is very low and amounts to 0.18 – 0.2 ml per tooth root. A tooth with 3 roots therefore requires a maximum of 0.6 ml.

The HSW Soft-Ject® dosing wheel syringe is compatible with all standardised glass cartridges with a volume of 1.7 ml or 1.8 ml. We recommend use of suitable short injection needles with short bevel, eg. HSW Fine-Ject® Dental 0.30 x 13mm / 30G x 1/2", or HSW Fine-Ject® r1Lax needle.

#### Preparation of device

- Use dosing wheel to insert toothed piston rod and move it back by turning the wheel into the instrument holder.
- Insert cartridge holder with instrument holder by locking the bayonet lock, the connection „clicks“ in.
- Move the dosing wheel by using thumb or fore finger towards the needle connector. The toothed piston rod moves into the cartridge holder. Move the piston rod until there is a light resistance. The tip of the toothed piston rod now touches the gasket of the cartridge. The front end of the cartridge is now pushed into the needle connector.
- Remove the back cap of the injection needle, insert the needle into the needle connector and screw the needle firmly onto the connector.
- Slowly move the dosing wheel towards the needle connector until the anaesthetics drops out of the needle.

#### Indication for intraligamental anaesthesia (ILA)

- Preparations of cavities and crowns.
- Endodontic treatment with the possibility of re-injections into the pulp area.
- Extractions of single teeth and osteotomies.
- Differential diagnosis of non-distinctive pulpitis pain.
- Treatment of patients with unstable bloodstream (due to the very low amount of anaesthetics required)
- Pain elimination for risk patients after suffering from heart attacks or cardiac bypasses
- Treatment of children or disabled patients for risk reduction of post-operative bite injuries. There are reports of enamel hypoplasia in permanent teeth induced by periodontal ligament anaesthesia of primary teeth.(Brännström, 1984).
- Patients suffering from haemorrhagic diathesis and undergoing anticoagulant therapy.
- Finalisation of partially failed nerve block anaesthesia.

#### Performing intraligamental anaesthesia (ILA) with the HSW Soft-Ject® dosing wheel syringe

Following the determination of the individual patient indication, we recommend the following procedure:

- **Determine each point of injection** for each root of the respective tooth.
- **Insert needle tip** alongside the tooth into **periodontal gap** 1 - max. 3 mm. The needle tip need to sit firmly in the sulcus and need to remain in this position during the injection.
- The free hand may give support for the needle introduction.
- **Move the dosing wheel** by using thumb or fore finger towards the needle tip; the toothed piston rod moves the gasket of the cartridge and anaesthetics is injected.  
Note: one „complete actuating path“ of the dosing wheel refers to moving the wheel by 5 scale marks.
- Overcome the sensible backpressure created by the tissue during injection.
- **Maintain the injection pressure** throughout the entire injection procedure.  
Each increment on the dosing wheel represents approximately injection of 0.02ml, depending on the manufacturing tolerances of the cartridge. (5 increments = 0.1 ml = 10 seconds)
- **Inject** approx. 0.18 – 0.2 ml of anaesthetics for each root of the tooth.
- **Injection time** requires a minimum of **≥ 20 seconds per root of the tooth**, for teeth with three roots > 60 seconds.
- Slow injection results in a gradual resorption of the anaesthetics into the spongy alveolar bone tissue.
- **Controlled and slow injection** avoids a **pressure pain** post analgesia
- **Anaesthesia** occurs instantly and is distractive after 30 seconds. After approx. 30 mins, the sensibility is fully re-established.

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