








**HENKE
SASS
WOLF**






Hysteroscopes 2.7 mm

Description		Part number	Working length	Total length
Hysteroscope 2.7 mm autoclavable 0°		7199202700	270 mm	340 mm
Hysteroscope 2.7 mm autoclavable 30°		7199202730	270 mm	340 mm
Obturator for diagnostic sheath		7019020035		
Continuous-flow diagnostic sheath with 2 rotatable stopcocks (outer diameter 5.5 mm)		7019020033	224 mm	
Continuous-flow operative sheath with 2 rotatable stopcocks (outer diameter 6.1 mm), working channel for flexible instruments diameter 2.1 mm		7019020045	190 mm	

Hysteroscopes 2.9 mm

Description		Part number	Working length	Total length
 Hysteroscope 2.9 mm autoclavable 0°		8300009034	300 mm	362 mm
 Hysteroscope 2.9 mm autoclavable 12°		8300009024	300 mm	362 mm
 Hysteroscope 2.9 mm autoclavable 30°		8300009029	300 mm	362 mm

Hysteroscopes 4.0 mm

Description		Part number	Working length	Total length
 Hysteroscope 4.0 mm autoclavable 0°		8300003244	300 mm	355 mm
 Hysteroscope 4.0 mm autoclavable 12°		8300003259	300 mm	355 mm
 Hysteroscope 4.0 mm autoclavable 30°		8300003261	300 mm	355 mm
 Hysteroscope 4.0 mm autoclavable 70°		8300003799	300 mm	355 mm