

# HSW ECO-MATIC® Drencher

12,5 ml / 30 ml / 70 ml

Automatic self-filling syringe for oral application  
Jeringa automática para aplicación oral

- Especially designed for sheep and cattle or multiple dosing of livestock, optionally available as injector
- Ergonomic and user-friendly design
- High-strength materials guarantee continuous use as well as long service life at low weight
- Dosing wheel for setting the required dosage
- Expansion of our renowned product line of HSW ECO-MATIC® syringes

- Desarrollado especialmente para ovejas y vacuno, facultativamente utilizable como inyector
- Diseño ergonómico y absolutamente fácil de usar
- El uso de materiales altamente resistentes y permanentes representa robustez y longevidad
- Gran resistencia a pesar de un peso ligero
- Rueda dosificadora para ajuste de la dosificación deseada
- Ampliación de la familia de productos de jeringas conocida HSW ECO-MATIC®




## Included in delivery

HSW ECO-MATIC® Drencher complete:  
1 x Syringe  
1 x Tubing  
1 x Drench nozzle 11 cm

## Alcance del suministro

HSW ECO-MATIC® Drencher completo:  
1 x Jeringa  
1 x Tubo  
1 x Adaptador de tubo 11 cm



Version  
Modelo

Vol/ml

Item no.  
No. de artículo

Drencher

12,5 ml (0,5 - 12,5 ml)	8300017257
30 ml ( 5 - 30 ml)	8300026525
70 ml (7,5 - 70 ml)	8300029836

Different nozzle attachments are available depending on the application needs.

Partes delanteras diferentes adaptables según aplicación.

