

HSW ECO-MATIC®

1 ml / 2 ml / 5 ml

Automatic self-filling syringe with threaded bottle adapter for oral application
Jeringa automática con adaptador de botella roscado para aplicación oral

- Compact bottle mount drencher
 - Durable, metal Luer Lock fitting
 - Accurate, continuously variable dose setting
 - Exceptional in price and performance
 - Ergonomic design for ease of use
 - Secure bottle fitting
 - Attaches to 28 mm threaded bottles
-
- Aplicador compacto para el montaje de una botella roscado
 - Cono Luer Lock de metal antiabrasión
 - Ajuste preciso y continuo de dosificación
 - Precio y prestaciones excepcionales
 - Diseño ergonómico, manejo muy sencillo
 - Fijación de botella segura
 - Se conecta a botellas roscadas de 28 mm



Nozzles are available separately, see accessories (p. 29)
Adaptador de tubo se vende por separado, véase accesorios (p. 29)

Included in delivery

HSW ECO-MATIC® complete:
1 x Syringe
1 x Threaded bottle adapter

Alcance del suministro

HSW ECO-MATIC® completo:
1 x Jeringa
1 x Adaptador de botella roscado

Version Modelo	Vol/ml	Item no. No. de artículo
Luer Lock	1 ml	8300018607
	2 ml	8300015213
	5 ml	8300018605



Threaded bottles available in different sizes: Botellas de rosca disponibles en colores y tamaños diferentes:

Version Modelo	Vol/ml	Item no. No. de artículo
Plastic bottle Botella de plástico	100 ml	8300019266
	250 ml	8300019267
	500 ml	8300017875

Lid sold separately /
El tapón se vende por separado: 8300017876

