



# ACCESSORIES

## FOR HSW ECO-MATIC® 6ml / 10ml

PRODUCT VIDEO:



### NOZZLES

The applicator has a standardized front part to which various nozzles can be attached (either for oral or pour-on application) as well as a Luer Lock attachment for injection needles.

This front part is used in the majority of our product range of HSW syringes allowing the different attachments to be interchangeable between (e.g. HSW ECO-MATIC® 6ml / 10ml, HSW ECO-MATIC® Drencher 12,5ml / 30ml, HSW DRENCH-MATIC®).

**a) Nozzle curved 5,5 inch (14 cm)**  
item code: 3600005605

**b) Nozzle curved 4,3 inch (11 cm)**  
item code: 3600003605

**c) Pour-On nozzle shower (plastic)**  
item code: 8300016387

**d) Nozzle Luer-Lock**  
item code: 8300016477

**e) Straight connecting tube, flat-jet (metal)**  
item code: 8300006381



**f) Straight connecting tube, flat-jet (plastic)**  
item code: 8300030096

**g) Nozzle curved (7,5 cm)**  
item code: 8300060721

### COLLARS FOR BACK BOTTLE MOUNT - VERSION



**h) Universal bottle collar white Ø 20mm**  
item code: 8300024134

**i) Universal bottle collar blue Ø 30mm**  
item code: 8300024136


**j) Universal bottle collar green Ø 33mm**  
item code: 8300024137

**k) Bottle collar for threaded bottles**  
item code: 8300016988

**l) Tubing collar 4-6mm**  
item code: 8300035585

## DRAW-OFF CAPS

- m) Draw-Off Cap white Ø 20 mm  
item code: 3600000020
- n) Draw-Off Cap blue Ø 30 mm  
item code: 3600000030
- o) Draw-Off Cap green Ø 33 mm  
item code: 3600000033



## HSW BOTTLE HOLDERS

- Arm model for 50 ml - 100 ml vaccine bottles  
item code: 8300027748
- Arm model for 250 ml vaccine bottles  
item code: 8300027749
- Arm model for 500 ml - 1000 ml vaccine bottles  
item code: 8300027747
  
- Belt / hanging model\* for  
50 ml - 100 ml vaccine bottles  
item code: 8300027746
- Belt / hanging model\* for  
250 ml vaccine bottles  
item code: 8300027750
- Belt / hanging model\* for  
500 ml – 1000 ml vaccine bottles  
item code: 8300027751



\*Belt/hanging model can be used on belt as shown or hanging around neck/shoulder.

## Following accessories are already included per item:

Item Number	Name	Included Accessories
8300059185	HSW ECO-MATIC® 6ml Injector Back Bottle Mount Universal	Nozzle Luer-Lock: 8300016477 Universal bottle collar Ø 20mm, Ø 30mm, Ø 33mm: 3600000020, 3600000030, 3600000033
8300059187	HSW ECO-MATIC® 10ml Injector Back Bottle Mount Universal	Nozzle Luer-Lock: 8300016477 Universal bottle collar Ø 20mm, Ø 30mm, Ø 33mm: 3600000020, 3600000030, 3600000033
8300059213	HSW ECO-MATIC® 6ml Drencher Tube Attachment	Nozzle curved 4,3 inch (11 cm): 3600003605
8300059215	HSW ECO-MATIC® 10ml Drencher Tube Attachment	Nozzle curved 4,3 inch (11 cm): 3600003605