



ACCESSORIES

FOR HSW ECO-MATIC® 6ml / 10ml



NOZZLES

The applicator has a standardized front part to which various nozzles can be attached (either for oral or pour-on application) as well as a Luer Lock attachment for injection needles.

This front part is used in the majority of our product range of HSW syringes allowing the different attachments to be interchangeable between (e.g. HSW ECO-MATIC® 6ml / 10ml, HSW ECO-MATIC® Drencher 12,5ml / 30ml, HSW DRENCH-MATIC®).

- a) **Nozzle curved 5,5 inch (14 cm)**
item code: 360005605
- b) **Nozzle curved 4,3 inch (11 cm)**
item code: 360003605
- c) **Pour-On nozzle shower (plastic)**
item code: 8300016387
- d) **Nozzle Luer-Lock**
item code: 3600030619
- e) **Straight connecting tube, flat-jet (metal)**
item code: 8300006381



- f) **Straight connecting tube, flat-jet (plastic)**
item code: 8300030096
- g) **Nozzle curved (7,5 cm)**
item code: 8300060734


COLLARS FOR BACK BOTTLE MOUNT - VERSION



- h) **Universal bottle collar white Ø 20mm**
item code: 8300024134
- i) **Universal bottle collar blue Ø 30mm**
item code: 8300024136
- j) **Universal bottle collar green Ø 33mm**
item code: 8300024137
- k) **Bottle collar for threaded bottles**
item code: 8300016988
- l) **Tubing collar 4-6mm**
item code: 8300035585

DRAW-OFF CAPS

- m) Draw-Off Cap white Ø 20 mm
item code: 360000020
- n) Draw-Off Cap blue Ø 30 mm
item code: 360000030
- o) Draw-Off Cap green Ø 33 mm
item code: 360000033



HSW BOTTLE HOLDERS

- Arm model for 50 ml - 100 ml vaccine bottles
item code: 8300027748
- Arm model for 250 ml vaccine bottles
item code: 8300027749
- Arm model for 500 ml - 1000 ml vaccine bottles
item code: 8300027747
- Belt / hanging model* for 50 ml - 100 ml vaccine bottles
item code: 8300027746
- Belt / hanging model* for 250 ml vaccine bottles
item code: 8300027750
- Belt / hanging model* for 500 ml – 1000 ml vaccine bottles
item code: 8300027751



*Belt/hanging model can be used on belt as shown or hanging around neck/shoulder.

Following accessories are already included per item:

Item Number	Name	Included Accessories
8300059185	HSW ECO-MATIC® 6ml Injector Back Bottle Mount Universal	Nozzle Luer-Lock: 3600030619 Universal bottle collar Ø 20mm, Ø 30mm, Ø 33mm: 8300024134, 8300024136, 8300024137
8300059187	HSW ECO-MATIC® 10ml Injector Back Bottle Mount Universal	Nozzle Luer-Lock: 3600030619 Universal bottle collar Ø 20mm, Ø 30mm, Ø 33mm: 8300024134, 8300024136, 8300024137
8300059213	HSW ECO-MATIC® 6ml Drencher Tube Attachment	Nozzle curved 4,3 inch (11 cm): 3600003605
8300059215	HSW ECO-MATIC® 10ml Drencher Tube Attachment	Nozzle curved 4,3 inch (11 cm): 3600003605