

NEW! HSW ECO-MATIC® 6ml / 10ml

AUTOMATIC SELF-FILLING SYRINGE FOR INJECTION, DRENCH- OR POUR-ON APPLICATION



FLEXIBILITY



Standardized swivel nut for easy exchange of needles, nozzles and hooks with other products from the HSW product range

Different adaptors for different shaped and sized bottles, as well as for tubes

USERFRIENDLY

Ergonomic and compact design

Easy valve and O-ring service

DOSE SETTING

Exact dosage adjustment via dosing wheel (0.5ml steps, from 0.5ml)

RESISTANT

Premium quality plastic syringe with durable plastic barrel

TWO SET-UPS

Back Bottle Mount or Tubing Version

ITEM CODES

8300059185	HSW ECO-MATIC® 6ml Injector Back Bottle Mount Universal
8300059187	HSW ECO-MATIC® 10ml Injector Back Bottle Mount Universal
8300059213	HSW ECO-MATIC® 6ml Drencher Tube Attachment
8300059215	HSW ECO-MATIC® 10ml Drencher Tube Attachment

The table only represents an excerpt of possible combinations in terms of dosing, adapters and front parts

For further information please contact:

Henke-Sass, Wolf GmbH
Keltenstraße 1
78532 Tuttlingen – Germany
Tel. +49 (0)7462 9466 - 0
E-Mail: info@henkesasswolf.de
www.henkesasswolf.de

PRODUCT VIDEO:

