



ACCESSORIES

FOR HSW ECO-MATIC® 6ml / 10ml

NOZZLES

The applicator has a standardized front part to which various nozzles can be attached (either for oral or pour-on application) as well as a Luer Lock attachment for injection needles.

This front part is used in the majority of our product range of HSW syringes allowing the different attachments to be interchanged between (e.g. HSW ECO-MATIC® 6ml / 10ml, HSW ECO-MATIC® Drencher 12,5ml / 30ml, HSW DRENCH-MATIC®).

a) Nozzle curved 5,5 inch (14 cm)
item code: 3600005605

b) Nozzle curved 4,3 inch (11 cm)
item code: 3600003605

c) Pour-On nozzle shower (plastic)
item code: 8300016387

d) Nozzle Luer-Lock
item code: 8300016477

e) Straight connecting tube, flat-jet (metal)
item code: 8300006381



f) Straight connecting tube, flat-jet (plastic)
item code: 8300030096

g) Nozzle curved (7,5 cm)
item code: 8300060721

COLLARS FOR BACK BOTTLE MOUNT - VERSION



h) Universal bottle collar white Ø 20mm
item code: 8300024134

i) Universal bottle collar blue Ø 30mm
item code: 8300024136


j) Universal bottle collar green Ø 33mm
item code: 8300024137

k) Bottle collar for threaded bottles
item code: 8300016988

l) Tubing collar Ø 4-6mm
item code: 8300035585

DRAW-OFF CAPS

- m) Draw-Off Cap white Ø 20 mm
item code: 3600000020
- n) Draw-Off Cap blue Ø 30 mm
item code: 3600000030
- o) Draw-Off Cap green Ø 33 mm
item code: 3600000033



HSW BOTTLE HOLDERS

- Arm model for 50 ml - 100 ml vaccine bottles
item code: 8300027748
- Arm model for 250 ml vaccine bottles
item code: 8300027749
- Arm model for 500 ml - 1000 ml vaccine bottles
item code: 8300027747
- Belt / hanging model* for
50 ml - 100 ml vaccine bottles
item code: 8300027746
- Belt / hanging model* for
250 ml vaccine bottles
item code: 8300027750
- Belt / hanging model* for
500 ml – 1000 ml vaccine bottles
item code: 8300027751



*Belt/hanging model can be used on belt as shown or hanging around neck/shoulder.

Following accessories are already included per item:

Item Number	Name	Included Accessories
8300059185	HSW ECO-MATIC® 6ml Injector Back Bottle Mount Universal	Nozzle Luer-Lock: 8300016477 Universal bottle collar Ø 20mm, Ø 30mm, Ø 33mm: 3600000020, 3600000030, 3600000033
8300059187	HSW ECO-MATIC® 10ml Injector Back Bottle Mount Universal	Nozzle Luer-Lock: 8300016477 Universal bottle collar Ø 20mm, Ø 30mm, Ø 33mm: 3600000020, 3600000030, 3600000033
8300059213	HSW ECO-MATIC® 6ml Drencher Tube Attachment	Nozzle curved 4,3 inch (11 cm): 3600003605
8300059215	HSW ECO-MATIC® 10ml Drencher Tube Attachment	Nozzle curved 4,3 inch (11 cm): 3600003605