

# HSW ECO-MATIC® 0,3 ml / 0,5 ml / 1 ml / 2 ml / 5 ml

Automatic self-filling syringe for injections with Bottle Mount Universal Collars  
 Jeringa automática para inyección con kit adaptador universal

- Color coded bottle collars for easy reference and identification of the right adapters for the different sized bottles
  - Easy use, durable design, safe and secured bottle fitting
  - For use with all medication bottles of different sizes and shapes (10, 50, 100, 200, 250 or 500 ml and larger) depending on the collar used: white Ø 20 mm, blue Ø 30 mm, green Ø 33 mm
- 
- Código de color de adaptadores de botella para identificación y referencia fácil
  - Uso fácil, diseño durable, colocación seguro de la botella
  - Para el uso con todas botellas de medicamento de tallas y formas diferentes (10, 50, 100, 200, 250 o 500 ml y más grandes) dependiente de los adaptadores usados: blanco Ø 20 mm, azul Ø 30 mm, verde Ø 33 mm



## Included in delivery

HSW ECO-MATIC® complete:  
 1 x Syringe  
 1 x bottle collar Ø 20 mm  
 1 x bottle collar Ø 30 mm  
 1 x bottle collar Ø 33 mm

## Alcance del suministro

HSW ECO-MATIC® completo:  
 1 x Jeringa  
 1 x Adaptador universal Ø 20 mm  
 1 x Adaptador universal Ø 30 mm  
 1 x Adaptador universal Ø 33 mm



## Version Modelo

## Vol/ml

## Item no. No. de artículo

Luer Lock

0,3 ml (0,1 ml - 0,3 ml)  
 0,5 ml (0,1 ml - 0,5 ml)  
 1 ml (0,1 ml - 1,0 ml)  
 2 ml (0,2 ml - 2,0 ml)  
 5 ml (0,5 ml - 5,0 ml)

8300026241  
 8300037027  
 8300026242  
 8300026243  
 8300026244

