

HSW MICRO-MATIC®

0,025 ml / 0,05 ml / 0,1 ml

Automatic self-filling syringe for micro-injections
 Jeringa automática para microinyecciones

- Exact dosage of smallest injection quantities (0.025 ml / 0.05 ml / 0.1 ml)
- Designed for single or dual injection via one needle
- Optimized for high hand comfort administering low dose injections
- High-strength materials
- Long service life and salt water resistant
- Mechanical counter 0 – 999.999 (automatic reset afterwards)
- Dosificación exacta de cantidades de inyecciones muy pequeñas (0,025 ml / 0,05 ml / 0,1 ml)
- Diseñado para inyección simple o doble a través de una aguja
- Optimizado para el confort de manos en administrando inyecciones de dosis bajas
- Materiales altamente resistentes
- Alta durabilidad, también en cuanto al contacto con agua salada
- Contador mecánico 0 – 999.999 (después reset automático)



Included in delivery (Double)

- 1 x Syringe incl. mechanical counter
- 2 x Silicone tubing 1,5 m
- 1 x Spare parts set
- 1 x Needle Guard

Alcance del suministro (Double)

- 1 x Jeringa con contador mecánico
- 2 x Tubo de silicona 1,5 m
- 1 x Kit de repuesto
- 1 x Needle guard

Included in delivery (Single)

- 1 x Syringe incl. mechanical counter
- 1 x Silicone tubing 1,5 m
- 1 x Spare parts set
- 1 x Needle Guard

Alcance del suministro (Single)

- 1 x Jeringa con contador mecánico
- 1 x Tubo de silicona 1,5 m
- 1 x Kit de repuesto
- 1 x Needle guard

		Item no. No. de artículo
HSW MICRO-MATIC® Double	0,05 / 0,025 ml	8300054350
HSW MICRO-MATIC® Double	0,05 / 0,05 ml	8300054349
HSW MICRO-MATIC® Double	0,1 / 0,05 ml	8300054347
HSW MICRO-MATIC® Double	0,1 / 0,1 ml	8300054348
HSW MICRO-MATIC® Single	0,05 ml	8300054345
HSW MICRO-MATIC® Single	0,1 ml	8300054346



Extra short needle sizes:

Tamaños de agujas extra cortas:

HSW-ECO®	0,70 x 3 mm = 22G x 1/8"
HSW-ECO®	0,70 x 4 mm = 22G x 1/8"
HSW-ECO®	0,60 x 3 mm = 23G x 1/8"
HSW-ECO®	0,60 x 4 mm = 23G x 1/8"
HSW-ECO®	0,60 x 5 mm = 23G x 1/5"

Cannulas on page 35 / Agujas en la página 35

