

# HSW MULTI-MATIC®

25 ml / 50 ml

Repeater syringe for multiple dosing of live stock  
 Jeringa dosificadora para vacunación en serie

- Durable pressure-cast aluminium handle
    - exceptional durability, yet lightweight
  - Durable medical grade plastic barrel
    - with metal needle fitting
    - biocompatible – autoclavable (134°C/273°F)
    - chemically resistant to all known medication
    - UV-resistant and break-proof
  - Exceptionally accurate dosing
    - tolerance < 1%
  - Easy dose setting
    - via 0,5ml (25 ml) or 1ml (50 ml) dosage wheel
  - Streamlined design
    - few spare parts required
  - Ergonomic design for ease of use
- 
- Mango de aluminio inyectado de alta resistencia
    - gran solidez con peso reducido
  - Cilindro de plástico de alta rigidez
    - con cono metálico de la aguja
    - biocompatible – apto para autoclave (134°C/273°F)
    - inalterable en contacto con todos los medicamentos conocidos
    - resistente a las radiaciones UV e irrompible
  - Gran exactitud en la dosificación
    - tolerancia < 1%
  - Fácil dosificación
    - mediante rueda dosificadora de 0,5 ml (25 ml) o de 1 ml (50 ml)
  - Fácil desmontaje para mantenimiento
    - pocas piezas de repuesto necesarias
  - Diseño ergonómico, manejo muy sencillo



**Included in delivery**

HSW MULTI-MATIC® complete:  
 1 x Syringe

**Alcance del suministro**

HSW MULTI-MATIC® completo:  
 1 x Jeringa

Version  
 Modelo

Vol/ml

Item no.  
 No. de artículo

Luer Lock

25 ml  
 50 ml

3625030051  
 3650030051