


HSW SAFE-MATIC®

2 ml

Automatic self-filling syringe for injections with safety features

Jeringa automática para inyección con características de seguridad

- 3 step safety mechanism to avoid accidental self-injections
 - Most suitable for medicines which could be harmful for the user
 - Syringe can be used universally with bottle adapter or tubing attachment
 - Medication bottles can be replaced quickly and easily (50 ml and 100 ml)
 - Durable and precise metal piston
 - Dosing fixed at 2 ml (not adjustable)
-
- 3 mecanismos de seguridad para evitar accidentes en la inyección
 - Más adecuado para los medicamentos que pueden ser nocivos para el usuario
 - La jeringa se puede utilizar con adaptador para botellas o con tubo de ascensión
 - Las botellas (50ml y 100 ml) se pueden cambiar fácilmente y rápidamente
 - Émbolo de metal preciso y duradero
 - Dosis fija 2 ml (sin graduación)



Included in delivery

HSW SAFE-MATIC® complete:
1 x Syringe
1 x Ascending tube 18 cm
1 x Tubing

Alcance del suministro

HSW SAFE-MATIC® completo:
1 x Jeringa
1 x Tubo de ascensión de 18 cm
1 x Tubo

Version

Modelo

Luer Lock

Vol/ml

2 ml

Item no.

No. de artículo

8300006508

The bottle adapter can be ordered separately.

El adaptador para botellas se puede suministrar por separado.