



# EPIG®

The first all-in-one handheld device for intramuscular needle-free vaccination of piglets

The EPIG® is a smart needle-free device for intramuscular vaccination of piglets, offering smooth and painless injections of 1 ml or 2 ml without recoil. With smart sensor feedback, Bluetooth connectivity and battery operation, it ensures high vaccination quality, improved animal welfare and reduced cross-contamination.

## ● ACCURATE AND VERSATILE DOSING

- highest dosing accuracy due to an integrated sensor control detecting and reporting air inside the liquid circuit
- dosage adjustable to 1 ml and 2 ml, vial sizes between 20 ml and 500 ml, compatible with most common intramuscular vaccines

## ● EFFICACY THROUGH SMART FEEDBACK AND EVALUATION

- optimum vaccination results thanks to direct feedback on each injection from multicolored LED bars on each side of the device
- integrated Bluetooth enables the evaluation of vaccination data after use (EPIG® App required, available for iOS and Android)

## ● SIMPLE AND USER-FRIENDLY

- no external power sources required - insert battery and go
- easy set-up, maintenance and cleaning procedures for long service life and highest biosecurity standards

## ● DURABLE AND SOPHISTICATED DESIGN

- robust medical grade components ensure splash water tightness (resistant acc. IP54)
- increasing vaccination quality and piglet welfare while reducing cross-contamination



### INTENDED USE:

Vaccination of piglets around weaning until the end of nursery

HSW ITEM CODE	ITEM DESCRIPTION
8300052120	EPIG® device (incl. two batteries and accessory set)