



HSW STATIONARY MOUNT FOR EPIG® DEVICE

FOR MORE INDEPENDENCE IN VACCINATION !

The stationary mount enables the independent performance of the vaccination process by lifting piglets up. For this the EPIG® is switched in stationary mode and attached to the stationary mount. The stationary mount can be placed securely and flexibly on any pen wall in the stable. To start the vaccination process a piglets neck has to be pushed against the EPIG® twice in quick succession.



YOUR BENEFITS

1) INDEPENDENCE

What normally requires two or more well-coordinated workers can be done by just one person with the help of the stationary mount. With this extension, not only more people can work at the same time, but also in their individual operating speed.

2) FLEXIBILITY

The mounting mechanism with quick release allows the operator to easily and quickly attach the stationary mount to any pig pen in the stable. This enables a flexible and unhindered movement of the workers.

3) INDIVIDUALITY

The height of the stationary mount is fully adaptable to the size of the operator and height of the pig pens. To ensure an ideal vaccination the angle of the head of the stationary mount, which holds the EPIG®, can be adjusted to the height of the device.

4) SAFETY

The EPIG® device contains a stationary mode which prevents unintended injections due to the pulled trigger. In this mode it is not possible to self-inject oneself by accidentally pushing the nozzle, because it has to be pushed twice in a short period of time to initiate the vaccination.



HSW ITEM CODE	ITEM DESCRIPTION
8300056421	Stationary mounting device for needle-free device (cpl.)