



# HSW Power-VACC

## Battery-Powered Vaccination Device

The HSW Power-VACC is a battery-powered injection device designed for intradermal, subcutaneous, intramuscular and nasal applications. It enables fast and precise injections for various animal species, ensuring reliable and accurate vaccination through its multifunctional application technology.

- **FAST, PRECISE AND VERSATILE**
  - for intradermal, subcutaneous, intramuscular and nasal administration
  - interchangeable cylinders: flexibility with 2 ml and 5 ml cylinders
  - fast and accurate injections for all animal species thanks to cutting-edge technology
- **BATTERY-POWERED**
  - high performance with battery life of up to 10,000 injections (depending on your settings)
  - compatible with standard lithium-ion batteries and chargers, which can be purchased locally
- **ERGONOMIC AND DURABLE DESIGN**
  - ergonomic handle ensures no fatigue, even during long periods of use
  - rugged design ideal for tough conditions, ensuring reliable performance in daily use
  - IP56 protection rating, easy to clean and optimal protection against dust and water
- **DIGITAL INTEGRATION**
  - HSW Connect App: detailed documentation, precise tracking and analysis of vaccination processes while keeping your data private and secure



PRODUCT VIDEO:



### CHOOSE BETWEEN...

HSW Power-VACC with intuitive control panel

HSW Power-VACC<sup>PRO</sup> with modern display and additional features,

e.g. visible counter, adjustable vial-level counter and position sensor to prevent air bubbles

HSW ITEM CODE	ITEM DESCRIPTION
8300064865	HSW Power-VACC Universal Version incl. battery (2 ml + 5 ml)
8300063994	HSW Power-VACC <sup>PRO</sup> Universal Version incl. battery (2 ml + 5 ml)
*MORE VERSIONS AVAILABLE ON REQUEST*	