








HSW Power-VACC / HSW Power-VACC^{Pro} Spare part price list 02/25



Item Code	Description	Picture
8300064001	Front part 2 ml, we recommend to exchange after every 25.000 injections	
8300064078	Front part 5 ml, we recommend to exchange after every 25.000 injections	
8300063576	Piston 2ml	
8300063583	Piston 5ml	
3602030080	O-Ring 2ml, we recommend to exchange after every 10.000 injections	0
3605030080	O-Ring 5ml, we recommend to exchange after every 10.000 injections	0
81402799	Basket	
8300056118	Snap bottle holder for spike	
8300034708	Bottle Mount Rubber Disc (standard)	
8300061637	Bottle Mount Rubber Disc XXL (for continuous use with threaded PET bottles)	
3600203033	Set Draw-Off-Caps \varnothing 20 + 30 + 33 mm (for connection of tubing with vial)	
8300016988	Threaded Bottle Adapter for Standard PET Threads	
8300019267	PET bottle 250ml (w/o cap)	
8300017876	PET bottle cap (white screw cap)	
8300064969	Set of three universal bottle collars (20 / 30 / 33 mm)	

8300035585	Tube connection for spike	
8300063836	Charger incl. cabel	
8300064460	Li-Ion-Akku 18650-2600mAh PCB protected	
8300063992	Battery cover	
8300036437	Silicone vent valve	
8300063958	pressure compensation element	