




HENKE
SASS
WOLF

HSW MULTI-MATIC® 50ml
Ersatzteilliste / Spare part list 02/25

Gültig ab / valid
01.03.2025



Nr.	Bezeichnung	Description	Appellation	HSW Cat. no.	MOQ
1	Plastikzylinder 50 ml LL	Plastic barrel 50 ml LL	Cylindre plastique 50 ml LL	3650030060	
2	O-Ring 50 ml	Washer 50 ml	Rondelle 50 ml	3650030080	10
3	Kolben+Kolbenstange 50 ml	Piston+piston rod 50 ml	Tige de piston+piston 50 ml	3650030070	
4	Rastblech komplett	Locking sleeve cpl.	Tôle d'enclench. cpl.	3625500220	10
5	Feder für Rastblech	Spring for locking sleeve	Ressort tôle d'enclench.	3625500210	10
6	Knopf schwarz	Knob black	Bouton noir	3625500100	10
7	Plastikzylinder 50 ml GW	Plastic barrel 50 ml thread	Cylindre plastique 50 ml à vis	3650050060	