



Nr.	Bezeichnung	Description	Appellation	No.	In set
1	Vorderteil LL	Nozzle Luer-Lock	Embout Luer-Lock	3000330809	
1b	Vorderteil Gewinde	Nozzle thread	Embout à vis	3000330819	
2	Dichtring	Washer	Rondelle	3000000090	1 pce.
3	Druckfeder	Compressing spring	Ressort à pression	3000000474	1 pce.
4	Metallkugel	Metal ball	Boule de métal	3000330001	2 pcs.
5	Zylinderbehälter	Barrel sleeve	Tube de cylindre	3001330804	
6	Führungshülse	Sleeve	Douille	8300016522	
7	Griff schwarz	Handle black	Poignée noire	3000000800	
8	Kolbenstange m. Griff	Piston rod w. handle	Tige de piston	8300004031	
9	Kolben	Piston	Piston	3001330805	
10	Einstellhülse	Adjusting sleeve	Manchon du disp. à dosage	8300016503	
11	Zylinderdichtung vorne	Washer	Rondelle	3001000814	1 pce.
12	Ventil-Druckfeder	Compressing spring	Ressort à pression	3000000085	1 pce.
13	Dichtring	Washer	Rondelle	3000000398	1 pce.
14	Druckfeder	Compressing spring	Ressort à pression	3000000813	
15	Rändelmutter	Knurled screw	Écrou moleté	8300016635	
16	Glaszylinder 1 ml	Glass barrel 1 ml	Cylindre en verre 1 ml	3001330820	
17	O-Ring	Washer	Rondelle	3000000607	1 pce.
18	Flaschenhalter ohne Korb	Bottle attachment w/o basket	Adaptateur bouteille sans panier	8300011133	
19	Korb für Flaschenhalter	Basket	Panier	81402799	
20	Ventilschraube	Valve screw	Vis tendeuse	8300016671	
21	Zylinderstift	Pin	Goujon	3000000812	
22	Zylinderdichtung hinten	Washer	Rondelle	3001000815	1 pce.
23a	Adapter f. Plastikflaschen	Adaptor plastic bottles	Adaptateur bouteille plastique	3000000461	
23b	Adapter f. Glasflaschen	Adaptor glass bottles	Adaptateur bouteille en verre	3000000471	
23	23a und 23b enthalten 24+25	23a and 23b contain 24 + 25	23a et 23b cont. 24 + 25		
24	Überwurfmutter	Fixing screw	Écrou	8300002045	
25	Schutzkappe	Closing cap	Capuchon à clôture	3000000403	
-	Ersatzteilset klein	Minor spare part kit	Set pièces de rechange min.	3003300100	
-	Ersatzteilset gross	Major spare part kit	Set pièces de rechange maj.	3003301101	
=	ET-Set klein+Zylinder	= minor kit+glass barrel	=set min.+cylindre en verre		