





**HENKE
SASS
WOLF**


Arthroscopes 1.9 mm

Description		Part number	Working length	Total length
Arthroscope 1.9 mm autoclavable 0°		7197201900	60 mm	110 mm
Arthroscope 1.9 mm autoclavable 30°		7197201930	60 mm	110 mm
Trocar sheath with 1 fixed stopcock, HSW-Lock (outer diameter 2.8 mm)		7012200160		
Trocar sheath with 2 fixed stopcocks 110°, HSW-Lock (outer diameter 2.8 mm)		7012200166		
Trocar triangular		7012200175		
Obturator blunt		7012200170		
Obturator conical blunt		7012200180		

Arthroscopes 2.7 mm

Description		Part number	Working length	Total length
Arthroscope 2.7 mm autoclavable 0°		7197202700	140 mm	190 mm
Arthroscope 2.7 mm autoclavable 30°		7197202730	140 mm	190 mm
Arthroscope 2.7 mm autoclavable 70°		7197202770	140 mm	190 mm
Trocar sheath with 2 rotatable stopcocks, HSW-Lock (outer diameter 4.0 mm)		7017200165		
Trocar sheath with 1 fixed stopcock, HSW-Lock (outer diameter 4.0 mm)		7017200160		
Trocar sheath with 1 rotatable stopcock, HSW-Lock (outer diameter 4.0 mm)		7017200150		
Trocar triangular		7017200175		
Obturator blunt		7017200170		
Obturator conical blunt		7017200180		

Arthroscopes 4.0 mm

Description		Part number	Working length	Total length
HD³ Arthroscope 4.0 mm autoclavable 0°		8300034737	185 mm	225 mm
HD³ Arthroscope 4.0 mm autoclavable 30°		8300034738	185 mm	225 mm
HD³ Arthroscope 4.0 mm autoclavable 70°		8300034739	185 mm	225 mm
Arthroscope 4.0 mm autoclavable 0°		8300009091	185 mm	225 mm
Arthroscope 4.0 mm autoclavable 30°		8300003242	185 mm	225 mm
Arthroscope 4.0 mm autoclavable 70°		8300009090	185 mm	225 mm
Quick-Lock trocar sheath with 2 rotatable stopcocks (outer diameter 6.0 mm)		8300037993		
Quick-Lock trocar sheath with 1 rotatable stopcock (outer diameter 6.0 mm)		8300037994		
Trocar triangular		8300037995		
Obturator blunt		8300037996		
Obturator conical blunt		8300037997		
Trocar sheath with 2 rotatable stopcocks, HSW-Lock (outer diameter 5.8 mm)		7017130165		
Trocar triangular		7017130175		
Obturator blunt		7017130170		
Obturator conical blunt		7017130180		