



Nr.	Bezeichnung	Description	Appellation	No.	MOQ
1	Überwurfmutter	Fixing screw	Écrou	8300017240	
2	Kunststoffzylinder 30 ml	Plastic barrel 30 ml	Cylindre en plastic 30 ml	3030000517	
3	Schlauch 1,7 m lang	Tube, length 1,7 m	Tuyau, 1,7 m	3561750600	
4	Ansatz Luer-Lock	Nozzle Luer-Lock	Embout Luer-Lock	3600030619	
5	Ansatzrohr gebogen 11 cm	Nozzle 11 cm	Embout 11 cm	3600003605	
-	Ansatzrohr gebogen 14 cm	Nozzle 14 cm	Embout 14 cm	3600005605	
6	Kolben 30 ml	Piston 30 ml	Piston 30 ml	8300040299	
7	Zylinderschraube	Screw for barrel	Écrou de cylindre	3000000531	
8	Feder f. Schlauch	Spring for tube	Ressort de traction our tuyau	82139999	5

**Ersatzteilsets / Spare part sets / Set pièces de rechange**

-	Ersatzteilset klein - bestehend aus Ventil vorne + Ventil hinten, O-Ring f. Kolben, O-Ring für Vorderteil, Kulissenstein, Filzring f. Kolben	consisting of front valve + back valve, O-ring for piston, O-ring for front part, sliding block, felt ring for piston	Set pièces de rechange petit - comprenant vanne distale + vanne proximale, joint torique pour piston, joint torique pour partie avant, coulisse, bague feutre pour piston	8300040463	
-	Ersatzteilset gross = Ersatzteilset klein + Zylinder	Major spare part set = Minor spare part set + cylinder	Set pièces de rechange grand = Set de pièces de rechange petit + cylindre	8300066922	
-	Ersatzteilset Zugfeder - bestehend aus oben gelb markierten Komponenten	Spare part set Tension Spring - consisting of components marked in yellow above	Set pièces de rechange pour ressort de traction - composé des éléments marqués en jaune ci-dessus	8300064638	