

Pneumatic Vaccinator - Spare parts / Consumables



Dripping chamber

Advised consumption: Once a day / Every vaccine exchange
Packaging: Single packed; sterilized

Item No.: 8300025431



HSW FINE-JECT® needle

Demand: min. one needle every 1.000 shots
Packaging: 100 pcs/box

Needle size: **0,90 x 25 mm / 20G x 1"**

Item No.: 4710009025



Spare part set

Demand: one pack every 100.000 shots

Item No. – **single**: 8300042994

Item No. – **double**: 8300042993



Tubing set

Item No.: 8300030544



Micro filter (air filter)

Demand: Once a month
Packaging: Single packed

Item No.: 8300049483



Compressed air coupler with hose adapter

Item No.: 8300031226



Cover for subcutaneous (neck) injection incl. trigger element

Item No.: 8300048428