


Eingeschränkte Verfügbarkeit / Supply not guaranteed


Nr.	Bezeichnung	Description	Appellation	No.	In set
1	Luer-Lock Ansatz	Nozzle Luer-Lock	Embout Luer-Lock	300000809	
2	Dichtring	Washer	Rondelle	300000090	1 pce.
3	Druckfeder	Compressing spring	Ressort à pression	300000474	1 pce.
4	Ventilkörper kpl. Viton	Valve with washer	Soupape avec rondelle	300000381	2 pcs.
5	Zylinderbehälter 2 ml	Barrel sleeve 2 ml	Tube de cylindre 2 ml	300200804	
6	Führungshülse	Sleeve	Douille	300200801	
7	Griff schwarz	Handle black	Poignée noire	300000800	
8	Kolbenstange	Piston rod	Tige de piston	300000818	
9	Kolben	Piston	Piston	300200805	
10	Einstellhülse	Adjusting sleeve	Manchon du disp. à dosage	8300016503	
11	Dichtring	Washer	Rondelle	300200814	1 pce.
12	Druckfeder	Compressing spring	Ressort à pression	300000085	1 pce.
13	Dichtring	Washer	Rondelle	3000000398	1 pce.
14	Druckfeder	Compressing spring	Ressort à pression	300000813	
15	Rändelmutter	Knurled screw	Écrou moleté	8300016635	
17	Glaszylinder 2 ml	Glass barrel 2 ml	Cylindre en verre 2 ml	300200820	
18	O-Ring	Washer	Rondelle	3000000607	1 pce.
19	Schlauchansatz	Tubing attachment	Raccord de tuyau	8300016643	
20	Griffstück kpl.	Crutch handle	Béquille	300000808	
21	Ventilschraube	Valve screw	Vis tendeuse	8300016670	
22	Zylinderstift	Pin	Goujon	300000812	
23	Dichtring	Washer	Rondelle	300200815	1 pce.
24	O-Ring	Washer	Rondelle	300200810	1 pce.
25	Steigrohr 18 cm	Ascending tube 18 cm	Tube à l'ascension 18 cm	3018760000	
26	Silikonschlauch 3 x 2	Silicone tube	Tuyau silicone	3532000600	
	Ersatzteilset 2 ml	Spare part set 2 ml	Set pièces de rechange 2 ml	300000200	
-	Vorderteil Gewinde	Nozzle thread	Embout à vis	300000819	