



Nr.	Bezeichnung	Description	Appellation	Art. No.	MOQ
1	Griff mit Kopf	Handle with head	Poignée avec tête	8300016539	
2	Griff	Handle	Poignée	8300016532	
3	Schlußschraube	Fixing screw	Vis de fixation à clôture	8300014832	5
4	Druckfeder	Compressing spring	Ressort à pression	3099000127	5
5	Bolzen	Pin	Goujon	3099000101	5
6	Druckfeder	Compressing spring	Ressort à pression	3099000108	5
7	Schaftschraube	Screw	Vis	3099001118	5
8	Knopf	Knob	Bouton	8300014831	5
9	Teilscheibe	Dosagesetting device	Disque diviseur	8300016664	
10	Fensterhülse	Metal tube	Tube en métal	8300016510	
11	Glaszylinder	Glass barrel	Cylindre	3010000070	5
12	Kolbenstange	Piston rod	Tige de piston	8300016604	
13	Zugstange komplett	Pullbar rod complete	Tige à tirer, cpl.	8300016687	
14	Kolbenscheibe	Piston plate	Disque de piston	3020000259	5
15	Vorderteil Luer-Lock	Nozzle Luer-Lock	Embout Luer-Lock	8300016675	
-	Vorderteil Gewinde	Nozzle Thread	Embout à vis	3010057500	

Ersatzteilsets / Spare Part Sets / Set pièces de rechange

-	Ersatzteilset 1, bestehend aus: Vorderer Dichtring, Kolben, hinterer Dichtring	Spare part set 1, consisting of: Washer rubber (front), Piston, Washer fibre (back)	Set pièces de rechange 1, composé de: Rondelle plat (avant), Rondelle (arrière)	3010000098	1,64
-	Ersatzteilset 2, bestehend aus: Klinkenfeder, Schaftschraube, Transporthaken	Spare part set 2, consisting of: Ratched spring, screw, locker arm ratchet	Set pièces de rechange 2, composé de: Ressort loquet, Vis, Crochet de transport	8300064079	4,12