



Nr.	Bezeichnung	Description	Appellation	HSW Cat. no.	In set
1	Griff mit Kopf	Handle with head	Poignée avec tête	8300016539	
2	Griff	Handle	Poignée	8300016532	
3	Schlußschraube	Fixing screw	Vis de fixation à clôture	8300014832	
4	Druckfeder	Compressing spring	Ressort à pression	3099000127	
5	Bolzen	Pin	Goujon	3099000101	
6	Druckfeder	Compressing spring	Ressort à pression	3099000108	
7	Schaftschraube	Screw	Vis	3099001118	
8	Knopf	Knob	Bouton	8300014831	
9	Teilscheibe	Dosagesetting device	Disque diviseur	8300016664	
9a	Sprengring	Spring-lock washer	Bague d'arrêt	85198199	
10	Dichtring	Washer fibre	Rondelle	3010000122	1 pce.
11	Fensterhülse	Metal tube	Tube en métal	8300016510	
12	Glaszylinder 10 ml	Glass barrel 10 ml	Cylindre en verre 10 ml	3010000070	
13	Kolbenstange	Piston rod	Tige de piston	8300016604	
14	Zugstange komplett	Pullbar rod complete	Tige à tirer, cpl.	8300016687	
15	Kolbenscheibe	Piston plate	Disque de piston	3020000259	
16	Kolben	Piston, silicone	Piston	3010000069	1 pce.
17	Dichtring	Washer rubber	Rondelle plat	3010000071	1 pce.
18	Vorderteil Luer-Lock	Nozzle Luer-Lock	Embout Luer-Lock	8300016675	
-	Vorderteil Recofix	Nozzle Recofix	Embout Recofix	3010017500	
-	Vorderteil Gewinde	Nozzle Thread	Embout à vis	3010057500	
19	Klinkenfeder	Ratchet spring	Ressort loquet	3099000112	
20	Schaftschraube	Screw	Vis	3099000107	
21	Transporthaken	Locker arm ratchet	Crochet de transport	3099000111	
-	Ersatzteilset	Spare part set	Set pièces de rechange	3010000098	