



Nr.	Bezeichnung	Description	Appellation	No.	In set	
					29	28
1	Vorderteil LL	Nozzle Luer-Lock	Embout Luer-Lock	3000330809		
1b	Vorderteil Gewinde	Nozzle thread	Embout à vis	3000330819		
2	Dichtring	Washer	Bague de garniture	3000000090	1	
3	Druckfeder	Compressing spring	Ressort à pression	3000000474	1	
4	Metallkugel	Metal ball	Boule de métal	3000330001	2	
5	Zylinderbehälter	Barrel sleeve	Tube de cylindre	8300008804		
6	Führungshülse	Sleeve	Douille	8300008938		1
8	Kolbenstange m. Griff	Piston rod w. handle	Barre de pist. avec poignée	8300008891		1
9	Kolben	Piston	Piston	8300009203		1
10	Druckfeder mit Nase	Compressing spring	Ressort à pression	8300009778		
11	Zylinderdichtung vorne	Washer in front	Joint cylindr. à l'avant	8300009208		
12	Ventil-Druckfeder	Compressing spring	Ressort à pression	3000000085	1	
13	Dichtring	Washer	Bague de garniture	3000000398	1	
14	Druckfeder	Compressing spring	Ressort à pression	3000000813		1
15	O-Ring	O-ring	Joint torique	85381399	1	1
16	Glaszylinder 2 ml	Glass barrel 2 ml	Cylindre en verre 2 ml	3002330820		
17	O-Ring	O-ring	Joint torique	3000000607	1	
18	Schlauchansatz	Tubing attachment	Raccord de tuyau	3000000806		
19	Korb	Basket	Panier	81402799		
20	Ventilschraube	Valve screw	Vis tendeuse	3000330458		
21	Hülse	Safety sleeve	Douille de sécurité	8300006509		
22	Kappe	Silicone cap	Capuchon de silicone	8300008822		
23	Steigrohr 18 cm	Ascending tube 18cm	Tube à l'ascension 18 cm	3018760000		
24	Silikonschlauch 3x2	Silicone tube	Tuyau silicone	3532000600		
25	Metallkugel	Metal ball	Boule de métal	3099000100		2
26	Druckfeder	Compressing spring	Ressort à pression	8300009542		
27	Flaschenhalter	Bottle attachment	Adaptateur bouteille	8300011133		
28	Ersatzteilset Griff	Spare part kit handle	Set pièces de rechange	8300013471		
29	Ersatzteilset klein*	Minor spare part kit*	Set pièces de rechange min.*	3003300200		
30	Ersatzteilset gross*	Major spare part kit*	Set pièces de rechange maj.*	3003301201		
	= ET-Set klein+Zylinder	= minor kit+glass barrel	=set min.+cylindre en verre			
	* = Spare part kit HSW UNI-MATIC®, which contains extra pieces not necessary for this syringe					