



| Nr. | Bezeichnung | Description | Appelation | Art. No. | MOQ |
|-----|--|---|---|------------|-----|
| 1 | Vorderteil LL | Nozzle Luer-Lock | Embout Luer-Lock | 3000330809 | |
| 2 | Führungshülse | Sleeve | Douille | 8300016522 | |
| 3 | Griff schwarz | Handle black | Poignée noire | 3000000800 | 5 |
| 4 | Kolbenstange m. Griff | Piston rod w. handle | Tige de piston complète | 8300004031 | |
| 5 | Kolben | Piston | Piston | 3001330805 | |
| 6 | Glaszylinder 1 ml | Glass barrel 1 cc | Cylindre en verre 1 ml | 3001330820 | |
| 7 | Flaschenhalter Universal | Bottle attachment w/o basket | Adaptateur bouteille sans panier | 8300038325 | |
| 8 | Zylinderstift | Pin | Goujon | 3000000812 | 5 |
| 9 | Set aus Universaladaptoren (Ø 20mm, 30mm, 33mm) | Set Universal Adapters (Ø 20mm, 30mm, 33mm) | Set d'adaptateurs universels (Ø 20mm, 30mm, 33mm) | 8300064969 | |
| 10 | Set aus Gewindeflaschenadapter, Schlauchadapter & Schlauch | Set of threaded bottle adapter, hose adapter & hose | Kit adaptateur pour bouteille filetée, adaptateur pour tuyau & tuyau en PVC | 8300064970 | |

Ersatzteilset / Spare Part Set / Set pièces de rechange

| | | | | | |
|---|---|---|---|------------|--|
| - | Ersatzteilset, bestehend aus: 2 x Ventilkugel, 2 x Ventulfeder, Dichtringen vorne und hinten f. Zylinder, Dichtring für LL Vorderteil | Spare parts set, consisting of: 2 x valve ball, 2 x valve spring, front and rear sealing rings for cylinder, sealing ring for LL front part | Set de pièces de rechange, comprenant : 2 x bille de soupape, 2 x ressort de soupape, bagues d'étanchéité avant et arrière pour cylindre, bague d'étanchéité pour partie avant LL | 3003300100 | |
|---|---|---|---|------------|--|