



Nr.	Bezeichnung	Description	Appellation	Art. No.	MOQ
1	Vorderteil LL	Nozzle Luer-Lock	Embout Luer-Lock	3000330809	
2	Führungshülse	Sleeve	Douille	8300016522	
3	Griff schwarz	Handle black	Poignée noire	3000000800	5
4	Kolbenstange m. Griff	Piston rod w. handle	Tige de piston complète	8300004031	
5	Kolben	Piston	Piston	3002330805	
6	Glaszylinder 2 ml	Glass barrel 2 cc	Cylindre en verre 2 ml	3002330820	
7	Flaschenhalter Universal	Bottle attachment w/o basket	Adaptateur bouteille sans panier	8300038325	
8	Zylinderstift	Pin	Goujon	3000000812	5
9	Set aus Universaladaptoren (Ø 20mm, 30mm, 33mm)	Set Universal Adapters (Ø 20mm, 30mm, 33mm)	Set d'adaptateurs universels (Ø 20mm, 30mm, 33mm)	8300064969	
10	Set aus Gewindeflaschenadapter, Schlauchadapter & Schlauch	Set of threaded bottle adapter, hose adapter & hose	Kit adaptateur pour bouteille fileté, adaptateur pour tuyau & tuyau en PVC	8300064970	

Ersatzteilsets / Spare Part Sets / Set pièces de rechange

-	Ersatzteilset, bestehend aus: 2 x Ventilkugel, 2 x Ventulfeder, Dichtringen vorne und hinten f. Zylinder, Dichtring für LL Vorderteil	Spare parts set, consisting of: 2 x valve ball, 2 x valve spring, front and rear sealing rings for cylinder, sealing ring for LL front part	Set de pièces de rechange, comprenant : 2 x bille de soupape, 2 x ressort de soupape, bagues d'étanchéité avant et arrière pour cylindre, bague d'étanchéité pour partie avant LL	3003300200	
---	---	---	---	------------	--