



Nr.	Bezeichnung	Description	Appellation	Art. No.	In set
1	Vorderteil LL	Nozzle Luer-Lock	Embout Luer-Lock	3000330809	
1b	Vorderteil Gewinde	Nozzle thread	Embout à vis	3000330819	
2	Dichtring	Washer	Rondelle	3000000090	1 pce.
3	Druckfeder	Compressing spring	Rondelle	3000000474	1 pce.
4	Metallkugel	Metal ball	Boule de métal	3000330001	2 pce.
5	Zylinderbehälter	Barrel sleeve	Tube de cylindre	8300003249	
6	Führungshülse	Sleeve	Douille	8300003176	
7	Griff schwarz	Handle black	Poignee noire	3000000800	
8	Kolbenstange m. Griff	Piston rod w. handle	Tige de piston complète	8300004031	
9	Kolben	Piston	Piston	8300003175	
10	Einstellhülse	Dosing sheath	Manchon du disp. à dosage	8300016503	
11	Zylinderdichtung vorne	Washer	Rondelle	8300003189	1 pce.
12	Ventil-Druckfeder	Compressing spring	Ressort à pression	3000000085	1 pce.
13	Dichtring	Washer	Rondelle	3000000398	1 pce.
14	Druckfeder	Compressing spring	Ressort à pression	3000000813	
15	Rändelmutter	Knurled screw	Écrou moleté	8300016635	
16	Glaszylinder 5 ml	Glass barrel 5 cc	Cylindre en verre 5 ml	8300005960	
17	O-Ring	Washer	Rondelle	3000000607	1 pce.
18	Schlauchansatz	Tubing attachment	Raccord de tuyau	8300016643	
20	Ventilschraube	Valve screw	Vis tendeuse	8300016671	
21	Zylinderstift	Pin	Goujon	8300003193	
22	Zylinderdichtung hinten	Washer	Rondelle	8300003190	1 pce.
23	Steigrohr 18 cm	Ascending tube 18 cm	Tube à l'ascension 18 cm	3018760000	
24	Silikonschlauch 3x2	Silicone tube	Tuyau silicon	3532000600	
	Ersatzteilset klein	Minor spare part kit	Set pieces de rechange min.	8300005958	
	Ersatzteilset gross	Major spare part kit	Set pieces de rechange maj.	8300011647	
	= ET-Set klein+Zylinder	= minor kit+glass barrel	= set min.+cylindre en verre		