



Nr.	Bezeichnung	Description	Appellation	Art. No.		In set
1	Vorderteil LL	Nozzle Luer-Lock	Embout Luer-Lock	3000330809		
1b	Vorderteil Gewinde	Nozzle thread	Embout à vis	3000330819		
2	Dichtring	Washer	Rondelle	3000000090		1 pce.
3	Druckfeder	Compressing spring	Ressort à pression	3000000474		1 pce.
4	Metallkugel	Metal ball	Boule de métal	3000330001		2 pcs.
5	Zylinderbehälter	Barrel sleeve	Tube de cylindre	3002330804		
6	Führungshülse	Sleeve	Douille	3002000801		
7	Griff schwarz	Handle black	Poignée noire	3000000800		
8	Kolbenstange m. Griff	Piston rod w. handle	Tige de piston complète	8300004031		
9	Kolben	Piston	Piston	3002330805		
10	Einstellhülse	Dosing sheath	Manchon du disp. à dosage	8300016503		
11	Zylinderdichtung vorne	Washer	Rondelle	3002000814		1 pce.
12	Ventil-Druckfeder	Compressing spring	Ressort à pression	3000000085		1 pce.
13	Dichtring	Washer	Rondelle	3000000398		1 pce.
14	Druckfeder	Compressing spring	Ressort à pression	3000000813		
15	Rändelmutter	Knurled screw	Écrou moleté	8300016635		
16	Glaszylinder 2 ml	Glass barrel 2 ml	Cylindre en verre 2 ml	3002330820		
17	O-Ring	Washer	Rondelle	3000000607		1 pce.
18	Schlauchansatz	Tubing attachment	Raccord de tuyau	8300016643		
20	Ventilschraube	Valve screw	Vis tendeuse	8300016671		
21	Zylinderstift	Pin	Goujon	3000000812		
22	Zylinderdichtung hinten	Washer	Rondelle	3002000815		1 pce.
23	Steigrohr 18 cm	Ascending tube 18 cm	Tube à l'ascension	3018760000		
24	Silikonschlauch 3x2	Silicone tube	Tuyau silicon	3532000600		
	Ersatzteilset klein	Minor spare part kit	Set pieces de rechange min.	3003300200		
	Ersatzteilset gross	Major spare part kit	Set pieces de rechange maj.	3003301201		
	= ET-Set klein+Zylinder	= minor kit+glass barrel	= set min.+cylindre en verre			