




HENKE
SASS
WOLF

HSW VET-MATIC® 10 ml
Ersatzteilliste / Spare part list 02/25

Gültig ab / valid
01.03.2025



Nr.	Bezeichnung	Description	Appellation	Art. No.	MOQ
1	Zylinderbehälter	Barrel sleeve	Tube de cylindre	8300016692	
2	Einschraubteil	Screw for front part	Écrou	3010000277	
3	Glaszylinder	Glass barrel	Cylindre en verre	3010000348	
4	Druckfeder	Compressing spring	Ressort à pression	3000000280	10
5	O-Ring	Washer	Rondelle	3010000373	10
6	Vorderteil Luer-Lock	Front part Luer-Lock	Embout Luer-Lock	3010030480	
7	Steigrohr 18 cm	Ascending tube 18 cm	Tube à l'ascension 18 cm	3018760000	
8	PVC-Schlauch	Rubber tube	Tuyau	3091820000	5

Ersatzteilsets / Spare part sets / Set pièces de rechange

9	Ersatzteilset vorderes Ventil	Spare part set front valve	Set pieces de rechange pour valve avant	8300064195	
10	Ersatzteilset hinteres Ventil	Spare part set back valve	Set pieces de rechange pour valve arrière	8300064197	