



Nr.	Bezeichnung	Description	Appellation	HSW Cat. no.	In set
1	Zylinderbehälter	Barrel sleeve	Tube de cylindre	8300016692	
2	Einschraubteil	Screw for front part	Écrou	3010000277	
3	Glaszylinder	Glass barrel	Cylindre en verre	3010000348	
4	Kolbenstange	Piston rod	Tige de piston	3000000275	
5	Druckfeder	Compressing spring	Ressort à pression	3000000280	
6	Einstellhülse	Adjusting sleeve	Manchon du disp. à dosage	8300016498	
7	Stellschraube	Fixing nut	Écrou de fixation	3010000378	
8	Schraubkappe	Screw cap	Capuchon de vis	3010000454	
9	Kolben	Piston	Piston	3010000479	
10	O-Ring	Washer	Rondelle	3010000373	2 pcs.
11	Ventil mit Dichtring	Valve with washer	Soupape avec rondelle	3000000281	2 pcs.
12	Ventil-Druckfeder	Valve spring	Ressort à pression	3000000474	1 pce.
13	Ventil-Druckfeder	Valve spring	Ressort à pression	3000000085	1 pce.
14	Vorderteil Luer-Lock	Front part Luer-Lock	Embout Luer-Lock	3010030480	
-	Vorderteil Gewinde	Front part thread	Embout à vis	3010050480	
15	O-Ring	Washer	Rondelle	3000000431	
16	Steigrohr 18 cm	Ascending tube 18 cm	Tube à l'ascension 18 cm	3018760000	
17	PVC-Schlauch	Rubber tube	Tuyau	3091820000	
-	Ersatzteilset	Spare part set	Set pièces de rechange	3010000099	



-	Ansatzrohr 11,4 cm	Nozzle 11,4 cm	Tuyau d'embase 11,4 cm	3600095605	
	"Mini-Drencher"	"Mini-Drencher"	"Mini-Administreur"		