



Nr.	Bezeichnung	Description	Appellation	No.	In set
1	Hülse	Sleeve	Tube en métal	8300016575	
2	Kunststoffzylinder 23 ml	Plastic barrel 23 ml	Cylindre en plastic 23 ml	3000000517	
3	Vorderteil	Front part	Tête de la seringue	8300016676	
4	Dichtring, rund	Washer, round	Rondelle, ronde	3000000526	
5	Mutter	Nut	Écrou de fixation	8300016630	
6	Ventil mit Dichtring, rund	Valve with washer	Soupape avec rondelle	3000000516	2 pce.
7	Ventil-Druckfeder	Valve spring	Ressort à pression	3000000622	1 pce.
8	Ansatzrohr 11cm	Nozzle 11cm	Embout 11cm	3000000605	
9	Überwurfmutter	Fixing screw	Écrou	8300016667	
10	Kolbenstange	Piston rod	Tige de piston	8300016615	
11	O-Ring	Washer	Rondelle	3000000606	1 pce.
12	Kolben	Piston	Piston	3000000780	
13	Griff, starr	Handle, rigid	Manche, fix	8300016542	
14	Stift	Splint	Goujon	3000000521	
15	Griff, beweglich	Handle, movable	Manche, mobile	8300016531	
16	Kulissenstein	Plastic insert	Glisseur	3600002532	
17	Achse	Pin	Axe	8300016472	
18	Einstellschraube	Fixation screw	Vis de calcage	8300016504	
19	Gegenmutter	Counter nut	Contre-écrou	8300016523	
20	Zugfeder	Draw spring	Ressort de traction	3000000616	
21	Spannschraube	Staining screw	Vis tendeuse	8300015858	
22	Spannbuchse	Locking bush	Douille tendeuse	3000001537	
23	Rändelmutter	Knurled screw	Écrou moleté	8300016633	
24	Zylinderschraube	Screw for barrel	Écrou de cylindre	3000000531	
25	Dichtring, rund	Washer, round	Rondelle, ronde	3000000607	
26	Schlauch 1,7 m lang	Tube, length 1,7 m	Tuyau, 1,7 m	3561750600	
27	Ventil-Druckfeder	Valve spring	Ressort à pression	3000000522	1 pce.
28	Scheibe	Disc	Disque	3000000539	
29	Sicherungsscheibe	Safety disc	Disque de protection	3000000621	
-	Ersatzteilset	Spare part set	Set pièces de rechange	3000000608	