

NEW! HSW F-Grip 0,5ml / 1ml / 2ml



AUTOMATIC SELF-FILLING SYRINGE FOR FATIGUE-FREE VACCINATIONS



FLEXIBILITY



Versatile use due to multiple models and accessories
bottle mount | tube feed | transformer

USERFRIENDLY

Ergonomic and compact design

Easy valve and O-ring service

DOSE SETTING

Exact dosage adjustment via dosing wheel

RESISTANT

Metal barrel body protects the cylinder

ECONOMICAL

Optimal price-benefit ratio and high-quality materials ensure competitive advantages

Available in three sizes and three set-ups

- 0,5ml / 1ml / 2ml
- Tubing, Bottle Mount or Transformer

Item no.	0,5ml	1ml	2ml
Tubing	8300032189	8300028679	8300028681
Bottle Mount	8300028710	8300028711	8300029815
Transformer	8300032142	8300032143	8300027782

For further information please contact:

Henke-Sass, Wolf GmbH
Keltenstraße 1
78532 Tuttlingen – Germany
Tel. +49 (0)7462 9466 - 0
E-Mail: info@henkesasswolf.de
www.henkesasswolf.de

