



FRESHLY HATCHED: OUR NEW SPRAY CABINET FOR DAY-OLD CHICKS

- Efficacy** energy-independent manual application for continuous operation
- Accurate** perfect droplet size due to an unique spray nozzle ensures even distribution of viral and coccidiosis vaccine
- Sophisticated** airflow in the 10-liter reservoir prevents sedimentation of coccidiosis vaccine
- User-friendly** self-explanatory operation as well as robust splash guard for user's safety
- Smart** easy set-up, maintenance and cleaning procedures for long service life and highest biosecurity standards
- Versatile** variable spray volume adjustment to either 16 or 22ml and suitable for crate depth of 400, 450 and 500mm
- Resistant** high-quality stainless steel components allow use with nearly every vaccine

Item No.:

8300048266 – HSW Spray Cabinet

For further information please contact:

Henke-Sass, Wolf GmbH
Keltenstraße 1
78532 Tuttlingen – Germany

Tel. +49 (0)7462 9466 380

E-Mail: michael.speck@henkesasswolf.de

www.henkesasswolf.de

