



THE FIRST ALL-IN-ONE HANDHELD NEEDLE-FREE DEVICE FOR INTRAMUSCULAR VACCINATION

- Efficacy** Optimum vaccination results thanks to direct feedback on each injection from multicolored LED bars on each side of the device
- Accurate** Highest dosing accuracy due to an integrated sensor control detecting and reporting air inside the liquid circuit
- Simple** No external power sources required - insert battery and go
- User-friendly** Easy set-up, maintenance and cleaning procedures for long service life and highest biosecurity standards
- Smart** Integrated Bluetooth enables the evaluation of vaccination data after use (EPIG® App required, available for iOS and Android)
- Resistant** Robust medical grade components ensure dust and water tightness (splash water resistant acc. IP54)
- Versatile** Dosage adjustable to 1ml and 2ml, vial sizes between 20 ml and 500 ml, compatible with most common intramuscular vaccines
- Sophisticated** Increasing vaccination quality and piglet welfare while reducing cross-contamination

Item No.:

8300052120 – EPIG® device (incl. two batteries and accessory set)

Intended use: Vaccination of piglets around weaning until the end of nursery

For further information please contact:

Henke-Sass, Wolf GmbH

Keltenstraße 1

78532 Tuttlingen – Germany

Tel. +49 (0)7462 9466 - 0

E-Mail: info@henkesasswolf.de

www.henkesasswolf.de

