



Dental Injection and Aspiration with Comfort

Our **HSW HENKE-DENT 500 with imperial (US) thread** featuring our **„Soft-Grip“** combines a classic cartridge syringe design with a new comfortable grip and comfort. Gripping plates and ring are coated with soft silicone that leads to a perfect handling of the

SOFT-GRIP blue (US):	8300042388
SOFT-GRIP red (US):	8300042387
SOFT-GRIP pink (US):	8300042385
SOFT-GRIP purple (US):	8300042386
SOFT-GRIP Petite blue (US):	8300042384
SOFT-GRIP Petite red (US):	8300042381
SOFT-GRIP Petite pink (US):	8300042383
SOFT-GRIP Petite purple (US):	8300042382



Product features:

- Reliable aspiration with harpoon aspiration device
- Comfortable handling through silicone coated gripping plates and ring
- Silicone available in different colours, OEM customers may choose from any colour to achieve optimal brand recognition
- Side-load for easy cartridge loading
- Fully autoclavable at 134 °C or 273 °F, extremely colour-fast

