

HSW ECO-MATIC® 0,3 ml / 0,5 ml / 1 ml / 2 ml / 5 ml

Automatic self-filling syringe for injections with Bottle Mount Basket

Jeringa automática para inyección con adaptador de botellas

- Premium-quality plastic syringe with durable plastic barrel
- Durable metal needle fitting
- Accurate, continuously variable dose setting
- Exceptional in price and performance
- Ergonomic, compact design for ease of use
- Durable plastic basket in two sizes, ensuring
 - secure attachment of medication bottle
 - bottle protected against shock or impact
 - medication bottles can be changed quickly and easily (50 ml to 250 ml)
- Jeringa de plástico de primera clase con cilindro de plástico de gran rigidez
- Cono de la aguja de metal antiabrasión
- Ajuste preciso y continuo de dosificación
- Precio y prestaciones excepcionales
- Diseño ergonómico, manejo muy sencillo
- Cesta en dos tallas plástico de alta resistencia, y por tanto
 - recepción segura de la botella del medicamento
 - protección de la botella frente a impactos y golpes
 - cambio fácil y rápido de las botellas de medicamento (50 a 250 ml)



Included in delivery

HSW ECO MATIC® complete:
1 x Syringe with bottle attachment

Alcance del suministro

HSW ECO-MATIC® completo:
1 x Jeringa con adaptador de botellas



Version

Modelo

Luer Lock

Vol/ml

0,3 ml (0,1 ml - 0,3 ml)
0,5 ml (0,1 ml - 0,5 ml)
1 ml (0,1 ml - 1,0 ml)
2 ml (0,2 ml - 2,0 ml)
5 ml (0,5 ml - 5,0 ml)

Item no.

No. de artículo

8300002189
8300034316
8300007181
3602130110
3605130110

