

HSW-PREMIUM®

Reusable veterinary needles

Agujas reutilizables para uso veterinario



- Distinctive bison bevel for long service life
- Tube made of spring-grade drawn chrome-nickel stainless steel
- Thick-walled needle tube with exceptional stability and elasticity
- Generously-shaped metal attachments
- Attachments are marked with mm and gauge information
- Direct removal from transparent container ensures reduced risk of injury and contamination
- Label features an ISO colour-code for quick identification of diameter
- Packaging: 10 items per case, 100 per box

- Ángulo pronunciado de biselaje para una mayor duración
- Tubo de acero inoxidable templado de cromo-níquel
- Tubo de la aguja de gran espesor con elevada estabilidad y elasticidad
- Amplia superficie de sujeción de metal
- Diámetro del tubo estampado en la superficie de sujeción, en mm y calibre
- Directamente extraíble de su contenedor transparente Bajo riesgo de lesión o contaminación
- Etiqueta provista de código de colores ISO para la rápida identificación del diámetro
- 10 unidades por envase, 100 por caja

Reusable needles (Luer Lock)

Agujas reutilizables (Luer Lock)

Dimensions	Colour-code	Item no.
Dimensiones	Código de color	No. de artículo
0,80 x 10 mm 21 G x 3/8"	10/08	3221080106*
0,90 x 10 mm 20 G x 3/8"	10/09	3221090106*
0,90 x 15 mm 20 G x 5/8"	15/09	3221090156*
1,00 x 10 mm 19 G x 3/8"	10/10	3221100106
1,00 x 15 mm 19 G x 5/8"	15/10	3221100156
1,00 x 20 mm 19 G x 3/4"	20/10	3221100206
1,00 x 25 mm 19 G x 1"	25/10	3221100256*
1,20 x 10 mm 18 G x 3/8"	10/12	3221120106*
1,20 x 15 mm 18 G x 5/8"	15/12	3221120156
1,20 x 20 mm 18 G x 3/4"	20/12	3221120206
1,20 x 25 mm 18 G x 1"	25/12	3221120256
1,20 x 35 mm 18 G x 1 3/8"	35/12	3221120356*
1,40 x 10 mm 17 G x 3/8"	10/14	3221140106*
1,40 x 15 mm 17 G x 5/8"	15/14	3221140156
1,40 x 20 mm 17 G x 3/4"	20/14	3221140206
1,40 x 25 mm 17 G x 1"	25/14	3221140256
1,40 x 30 mm 17 G x 1 1/4"	30/14	3221140306*
1,40 x 35 mm 17 G x 1 3/8"	35/14	3221140356
1,60 x 10 mm 16 G x 3/8"	10/16	3221160106*
1,60 x 15 mm 16 G x 5/8"	15/16	3221160156
1,60 x 20 mm 16 G x 3/4"	20/16	3221160206
1,60 x 25 mm 16 G x 1"	25/16	3221160256
1,60 x 30 mm 16 G x 1 1/4"	30/16	3221160306*
1,60 x 35 mm 16 G x 1 3/8"	35/16	3221160356
1,60 x 100 mm 16 G x 4"	100/16	3221161006**
1,80 x 15 mm 15 G x 5/8"	15/18	3221180156
1,80 x 20 mm 15 G x 3/4"	20/18	3221180206
1,80 x 25 mm 15 G x 1"	25/18	3221180256
1,80 x 30 mm 15 G x 1 1/4"	30/18	3221180306*
1,80 x 35 mm 15 G x 1 3/8"	35/18	3221180356
1,80 x 40 mm 15 G x 1 1/2"	40/18	3221180406*
2,00 x 15 mm 14 G x 5/8"	15/20	3221200156*
2,00 x 20 mm 14 G x 3/4"	20/20	3221200206
2,00 x 25 mm 14 G x 1"	25/20	3221200256
2,00 x 30 mm 14 G x 1 1/4"	30/20	3221200306*
2,00 x 35 mm 14 G x 1 3/8"	35/20	3221200356

Dimensions	Colour-code	Item no.
Dimensiones	Código de color	No. de artículo
2,00 x 40 mm 14 G x 1 1/2"	40/20	3221200406*
2,00 x 50 mm 14 G x 2"	50/20	3221200506**
2,20 x 25 mm -	25/22	3221220256*
2,20 x 30 mm -	30/22	3221220306*
2,20 x 35 mm -	35/22	3221220356*
2,40 x 25 mm 13G x 1"	25/24	3221240256*
2,40 x 30 mm 13G x 1 1/4"	30/24	3221240306*

Reusable needles (thread)

Agujas reutilizables (cierre roscado)

Dimensions	Colour-code	Item no.
Dimensiones	Código de color	No. de artículo
1,00 x 10 mm 19 G x 3/8"	10/10	3225100106
1,00 x 15 mm 19 G x 5/8"	15/10	3225100156
1,00 x 20 mm 19 G x 3/4"	20/10	3225100206
1,20 x 15 mm 18 G x 5/8"	15/12	3225120156
1,20 x 20 mm 18 G x 3/4"	20/12	3225120206
1,20 x 25 mm 18 G x 1"	25/12	3225120256
1,40 x 15 mm 17 G x 5/8"	15/14	3225140156
1,40 x 20 mm 17 G x 3/4"	20/14	3225140206
1,40 x 25 mm 17 G x 1"	25/14	3225140256
1,40 x 35 mm 17 G x 1 3/8"	35/14	3225140356
1,60 x 15 mm 16 G x 5/8"	15/16	3225160156
1,60 x 20 mm 16 G x 3/4"	20/16	3225160206
1,60 x 25 mm 16 G x 1"	25/16	3225160256
1,60 x 35 mm 16 G x 1 3/8"	35/16	3225160356
1,80 x 15 mm 15 G x 5/8"	15/18	3225180156
1,80 x 20 mm 15 G x 3/4"	20/18	3225180206
1,80 x 25 mm 15 G x 1"	25/18	3225180256
1,80 x 35 mm 15 G x 1 3/8"	35/18	3225180356
2,00 x 20 mm 14 G x 3/4"	20/20	3225200206
2,00 x 25 mm 14 G x 1"	25/20	3225200256
2,00 x 35 mm 14 G x 1 3/8"	35/20	3225200356
2,20 x 35 mm -	35/22	8300001897

