

HSW HENKE-JECT® Insulin syringes

0,5 ml / 1 ml

HSW HENKE-JECT® Jeringas de insulina

- Crystal clear cylinder with easy to read scale
- Piston, piston-rod and cylinder consisting of Polypropylene
- With integrated needle
- Free of Pyrogene and PVC
- Free of Latex
- Non-toxic
- Certified acc. to ISO 9626:1991/Amd 1:2001 and ISO 8537:2007
- Double o-ring for additional security
- Cilindro cristalino y escala bien visible
- Pistón y cilindro de polipropileno, tapón de pistón de caucho-poliisopreno sintético
- Con aguja integrada
- Sin pirógena y PVC
- Sin látex
- No tóxico
- Certificado según ISO 9626:1991/Amd 1:2001 y ISO 8537:2007
- Aro de guarnición doble para más seguridad



Version Modèle	Vol/ml	Item no. No d'article
0,5ml INS 40IU 29Gx1/2"	0,5 ml	8300034658
1ml INS 40IU 29Gx1/2"	1 ml	8300034660



Syringes of size 0.5 ml are delivered in a sterile polybag, 10 pcs. each.
Syringes of size 1 ml are delivered in a blister packaging.

Jeringas en tamaño 0,5 ml se suministrarán en bolso polietileno esteril (10 piezas).
Jeringas en tamaño 1 ml se suministrarán en envase blíster.

