



Nr.	Bezeichnung	Description	Appellation	No.	In set	
1	Hülse	Sleeve	Tube en métal	8300016576		
2	Kunststoffzylinder 30 ml	Plastic barrel 30 ml	Cylindre en plastic 30 ml	3030000517		
3	Vorderteil	Front part	Tête de la seringue	8300016678		
4	Dichtring, rund	Washer, round	Rondelle, ronde	3030000526		
5	Mutter	Nut	Écrou de fixation	8300016630		
6	Ventil mit Dichtring, rund	Valve with washer	Soupape avec rondelle	3000000516	2 pce.	
7	Ventil-Druckfeder	Valve spring	Ressort à pression	3000000622	1 pce.	
8	Ansatzrohr 11cm	Nozzle 11cm	Embout 11cm	3000000605		
9	Überwurfmutter	Fixing screw	Écrou	8300016667		
10	Kolbenstange 30 ml	Piston rod 30 ml	Tige de piston 30 ml	8300016616		
11	O-Ring	Washer	Rondelle	3030000606	1 pce.	
12	Kolben	Piston	Piston	3030000780		
13	Griff, starr	Handle, rigid	Mance, fix	8300016551		
14	Stift	Splint	Goujon	3000000521		
15	Griff, beweglich	Handle, movable	Manche, mobile	8300016531		
16	Kulissenstein	Plastic insert	Glisseur	3600002532		
17	Achse	Pin	Axe	8300016472		
18	Einstellschraube	Fixation screw	Vis de calcage	8300016504		
19	Gegenmutter	Counter nut	Contre-écrou	8300016523		
20	Zugfeder	Draw spring	Ressort de traction	3000000616		
21	Spannschraube	Staining screw	Vis tendeuse	8300015858		
22	Spannbuchse	Locking bush	Douille tendeuse	3000001537		
23	Rändelmutter	Knurled screw	Écrou moleté	8300016633		
24	Zylinderschraube	Screw for barrel	Écrou de cylindre	3000000531		
25	Dichtring, rund	Washer, round	Rondelle, ronde	3000000607		
26	Schlauch 1,7 m lang	Tube, length 1,7 m	Tuyau, 1,7 m	3561750600		
27	Ventil-Druckfeder	Valve spring	Ressort à pression	3000000522	1 pce.	
28	Scheibe	Disc	Disque	3000000539		
29	Sicherungs-scheibe	Safety disc	Disque de protection	3000000621		
-	Ersatzteilset	Spare part set	Set pièces de rechange	3030000608	2,12	