



Nr.	Bezeichnung	Description	Appellation	No.	In set
1	Überwurfmutter	Fixing screw	Écrou	8300017240	
2	Plastikzylinder 10 ml	Plastic barrel 10 ml	Cylindre plastique 10 ml	3610000517	
3	Überwurfmutter Griff	Fixing screw (handle)	Écrou (manche)	8300016668	
4	Kulissenstein	Plastic insert	Glisseur	3600002532	1 pce.
5	Schlauchanschluss	Tubing attachment	Raccord de tuyeau	8300040335	
6	PVC-Schlauch 1,7 m lang	PVC Tube length 1,7 m	Tuyau PVC 1,7 m	3561750600	
7	Zugfeder	Draw spring	Ressort de traction	3000000616	
8	Spannschraube	Staining screw	Vis tendeuse	8300015858	
9	Spannbuchse	Clamping bushing	Douille tendeuse	8300016656	
10	Scheibe	Disc	Disque	3000000539	
11	Rändelmutter	Knurled screw	Écrou moleté	8300016633	
12	Einstellschraube	Fixation screw	Vis de calcage	3600002534	
13	Rändelmutter	Knurled screw	Écrou moleté	8300016634	
14	Ansatz Luer-Lock	Nozzle Luer-Lock	Embout Luer-Lock	3600030619	
15	Ansatzrohr gebogen 11 cm	Drench nozzle 3,5"	Embout courbé 11 cm	3600003605	
-	Ansatzrohr gebogen 14 cm	Drench nozzle 5,5"	Embout courbé 14 cm	3600005605	
16	Kolbenstange 10 ml kpl.	Piston rod 10 ml cpl.	Tige de piston cpl. 10 ml	8300041497	
17	Griff 10 ml kpl.	Handle 10 ml cpl.	Manche fix cpl. 10 ml	8300041515	
18	O-Ring Kolben	Piston o-ring	Rondelle piston	3605030080	2 pcs.
19	O-Ring	O-ring housing rear valve	Rondelle	3600002606	1 pce.
20	O-Ring	O-ring LL+Drench nozzle	Rondelle	3000000519	1 pce.
21	O-Ring	O-ring mounting hub	Rondelle	3610002608	1 pce.
22	O-Ring	O-Ring	O-Ring	8300027215	1 pce.
23	Ventil vorne	Front valve	Soupape tête de la seringue	8300032159	1 pce.
24	Ventil hinten	Rear valve	Soupape arrière	8300032110	1 pce.
25	Feder f. Schlauch	Spring for tube	Ressort de traction pour tuyau	82139999	
26	O-Ring	O-Ring	O-Ring	3610002526	1 pce.
-	Ersatzteilset klein	Minor spare part set	Set pièces de rechange mini	8300040467	
-	Ersatzteilset gross	Major spare part set	Set pièces de rechange grand	8300040474	
=	ET-Set klein+Zylinder	=Minor set+plastic barrel	= set mini+cylindre plastique		