



Nr.	Bezeichnung	Description	Appellation	No.	In set
1	Hülse	Sleeve	Tube en métal	8300016576	
2	Kunststoffzylinder 30 ml	Plastic barrel 30 ml	Cylindre en plastic 30 ml	3030000517	
3	Vorderteil	Front part	Tête de la seringue	8300040294	
4	O-Ring	Washer	Rondelle	3030000526	
5	Mutter	Nut	Écrou de fixation	8300040295	
6	Ventil vorne kpl.	Valve cpl.	Soupape cpl.	8300032159	1 pce.
7	Ventil hinten kpl.	Valve cpl.	Soupape cpl.	8300042009	1 pce.
8	Ansatzrohr gebogen 11 cm	Nozzle 11 cm	Embout 11 cm	3600003605	
-	Ansatzrohr gebogen 14 cm	Nozzle 14 cm	Embout 14 cm	3600005605	
9	Überwurfmutter	Fixing screw	Écrou	8300017240	
10	Kolbenstange 30 ml	Piston rod 30 ml	Tige de piston 30 ml	8300016616	
11	O-Ring	Washer	Rondelle	8300027065	1 pce.
12	Kolben 30 ml	Piston 30 ml	Piston 30 ml	8300040299	
13	Griff, starr	Handle, rigid	Mance, fix	8300016551	
14	Stift	Splint	Goujon	3000000521	
15	Griff, beweglich	Handle, movable	Manche, mobile	8300016531	
16	Kulissenstein	Plastic insert	Glisseur	3600002532	
17	Achse	Pin	Axe	8300016472	
18	Einstellschraube	Fixation screw	Vis de calcage	8300016504	
19	Gegenmutter	Counter nut	Contre-écrou	8300016523	
20	Zugfeder	Draw spring	Ressort de traction	3000000616	
21	Spannschraube	Staining screw	Vis tendeuse	8300015858	
22	Spannbuchse	Locking bush	Douille tendeuse	3000001537	
23	Rändelmutter	Knurled screw	Écrou moleté	8300016633	
24	Zylinderschraube	Screw for barrel	Écrou de cylindre	3000000531	
25	O-Ring	Washer	Rondelle	3000000519	
26	Schlauch 1,7 m lang	Tube, length 1,7 m	Tuyau, 1,7 m	3561750600	
27	Scheibe	Disc	Disque	3000000539	
28	Sicherungsscheibe	Safety disc	Disque de protection	3000000621	
29	Filzring 30ml	Felt washer 30 ml	Bague feutre 30 ml	8300012205	1 pce.
30	O-Ring	Washer	Rondelle	8300027215	
31	Feder f. Schlauch	Spring for tube	Ressort de traction pour tuyau	82139999	
-	Ersatzteilset	Spare part set	Set pièces de rechange	8300040463	