

# ► HENKE-ASPIRATIONSSPRITZE

## nach Cook-Waite (CW Type)

1,8ml

*Zylinder-Ampullen-Spritze mit Aspiration*  
*Dental cartridge syringe with aspiration*



- Spritzendesign nach Cook-Waite
- seitliche Beladung
- benutzbar mit allen marktüblichen Zylinderampullen mit Vollstopfen
- Handstück und Ziehring frei drehbar für höchsten Anwenderkomfort
- Vorderteil verklebt und nicht abnehmbar
- zuverlässige Aspiration mit Pfeilspitze
- sichere Handhabung und leichte Reinigung
- autoklavierbar bei 134 °C oder 273 °F
  
- syringe design according to Cook-Waite
- side loading
- suitable for all commercially available dental cartridges with solid plug
- finger plates and thumb ring are completely rotatable for highest user comfort and flexibility
- front part permanently attached
- reliable aspiration with harpoon aspiration device
- secure handling and easy cleaning
- autoclavable at 134 °C or 273 °F



► Ausführung/Version	Vol./ml	Artikel-Nr./Item no.
Metrisches Gewinde/ metric thread	1,8 ml	2931018101
US-Gewinde/ imperial (US) thread	1,8 ml	2930718101

