

# ► HENKE-DENT 500 SOFT-GRIP®

1,8ml

## Zylinder-Ampullen-Spritze mit Aspiration Dental cartridge syringe with aspiration

- zur angenehmeren Handhabung Fingerstütze und Ziehring mit Silikon umspritzt
- verfügbar in 2 Größen (Standard und Petite) und in 4 Farben
- die Größe Petite zeichnet sich durch einen kleineren ovalen Ziehring und eine verkürzte Kolbenstange aus
- benutzbar mit allen marktüblichen Zylinderampullen mit Vollstopfen
- zuverlässige Aspiration mit Pfeilspitze
- seitliche Beladung
- Vorderteil verklebt und nicht abnehmbar
- sichere Handhabung und leichte Reinigung
- autoklavierbar bei 134°C oder 273°F
- the thumb ring and finger plates are silicone-covered to improve ergonomic handling and user comfort
- available in two sizes (standard and petite) and in 4 colors
- size petite is equipped with a smaller thumb ring and a reduced piston rod
- can be used with all commercially available cylinder ampoules with solid plug
- reliable aspiration with harpoon aspiration device
- side loading
- front part permanently attached
- secure handling and easy cleaning
- autoclavable at 134°C or 273°F

### ► **Lieferumfang**

1x Spritze

### ► **Included in delivery**

1x syringe

| ► <b>Ausführung/Version</b> | <b>Vol./ml</b> | <b>Artikel-Nr./Item no.</b> |
|-----------------------------|----------------|-----------------------------|
| Standard, blau/blue         | 1,8ml          | 8300010284                  |
| Standard, rot/red           | 1,8ml          | 8300010285                  |
| Standard, pink/pink         | 1,8ml          | 8300011934                  |
| Standard, violett/purple    | 1,8ml          | 8300011933                  |
| Petite, blau/blue           | 1,8ml          | 8300010276                  |
| Petite, rot/red             | 1,8ml          | 8300011936                  |
| Petite, pink/pink           | 1,8ml          | 8300010283                  |
| Petite, violett/purple      | 1,8ml          | 8300011935                  |

